

MAXI REST BARIATRIC BED

MAXI REST 4842 SERIES

Features

- Robust steel two-piece split-frame construction for ease of delivery and transport
- Bed deck width is expandable to 39", 42", or 48"
- Bed deck length is adjustable to 80" or 88"
- 5" Deluxe Caster Set includes 2 Total-Lock and 2 Steer-Lock for ease of maneuverability
- Standard 3-Function Positioning hand control
- Optional 5-Function Positioning hand control (Trendelenburg / Reverse Trendelenburg and Chair Positioning)
- Synchronized Hi/Lo Motors ensure bed remains level
- Battery backup standard for dependability
- Any height rolling position for ease of use
- Optional Intuitive side rail system for both head and foot sections
- Optional laminate wood grain or plastic molded head and foot board options
- Optional factory-installed weight scale system
- 800 lb (362.9 kg)¹ weight capacity¹ –EVENLY DISTRIBUTED

4842 Maxi Rest Bed Product Info



Options:

	Maxi Rest Bed Only	Split Frame	Rail Mounting at head and foot	Scale
GS4842D	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
GS4842DNS	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
GS4842DW	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

★ options available at additional cost

Specifications:

Series	Bed Deck Width Positions						Bed Deck Length Positions				Bed Deck Height Range				Weight ¹ Capacity EVENLY DISTRIBUTED		Standard	Optional	Motors	emergency Battery Backup
	inch		cm		inch		cm		inch		cm		inch		cm		lb	kg ¹	functions	features
	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum	minimum	maximum				
4842	39.00	99.06	42.00	106.68	48.00	121.92	80.00	203.20	88.00	223.52	13.50	34.29	27.50	69.85	800	362.87	3	5	4	<input checked="" type="checkbox"/>

★ options available at additional cost

† standard rounding method used

¹ maximum weight capacity - evenly distributed - refer to user manual

² accessories sold separately



MAXI REST 4842 SERIES



Maxi Rest Options and Accessories

sold separately at additional cost

- GL928000 38 inch (96.52 cm) Laminate Wood Grain Head and Foot Board Set
- GL929000 39 inch (99.06 cm) HDPE Plastic Head and Foot Board Set
- GL930000 48 inch (121.92 cm) HDPE Plastic Head and Foot Board Set
- GL927100 Deluxe Plastic Rail Set (LH/RH)
- GL927000 Standard Painted Head End Rails
- GL927020 Standard Painted Foot End Rails w/ 39, 42 & 48 inch Pull-outs
- GL93100 80 inch (203.20 cm) Frame Mount Full Length Trapeze [FITS 80 INCH ONLY]
- GL93180 88 inch (223.52 cm) Frame Mount Full Length Trapeze [FITS 88 INCH ONLY]
- GL93200 5-Function Hand Control (Trendelenburg/Rev Trendelenburg/Chair Positioning)
- GL93300 Full Width Wall Bumper
- GL93600L Egress (Assist) Bar and Left Hand Mounting Pullout
- GL93600R Egress (Assist) Bar and Right Hand Mounting Pullout
- GL94000 8 inch (20.32 cm) Bolster Mattress and Extension Pan Kit

GL929000

39 INCH HDPE PLASTIC HEAD AND FOOT BOARD SET



GL927000

STANDARD PAINTED RAIL SET SHOWN WITH HAND CONTROL



GL927100

DELUXE PLASTIC RAIL SET SHOWN WITH HAND CONTROL



Warranty

Main Frame and welds:..... Five years
Control Box and Actuator Motors:..... Two years
Hand Control Pendant, Staff Control, and Cabling:..... Two years
Head Board and Foot Board:..... One year
All other durable components not listed above:..... Two years


GL93200

5-FUNCTION TREND/REV
TREND HAND CONTROL



GL93100

FRAME MOUNT FULL LENGTH TRAPEZE
80 inch (203.20 cm) version shown



GS4824DW

BED WITH FACTORY INSTALLED SCALE SYSTEM



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.

graham
field