

graham field™

Weighing In on Graham-Field Offerings

Providing optimal care for the bariatric population presents a challenge for caregivers today. Statistics show that around 34% of the U.S. population is considered to be clinically obese. This epidemic has put a tremendous stress on the

health care system. Turning a patient over in bed, using a bedside commode, lateral transfers, and lifting someone from a bed to a chair are the everyday tasks associated with the patient care of an obese individual and can sometimes

require a specially designed product to help. Graham-Field offers a broad range of products designed with the larger individual in mind to help you focus more on the patient without the worries of a product to best suit their needs.

B700 BARIATRIC-BED

Item ABL-B700

A heavy-duty full-electric bed designed to support up to 600 pounds (patient weight). Its high-strength, split-frame make it durable and easy to deliver and transport.

Other features include:

- Grid sleep surface – ideal for heavier users to prevent bottoming out
- Overall 42" wide sleep surface (46.5" with rails)
- Powerful DC motors for head/foot adjustment as well as bed height

- Adjustment range from a low 17" to a high 27"
- Comes standard with emergency crank for backup
- Available with rails for enhanced patient safety
- Meets requirements for HCPCS E0301/E0303





**Bariatric Foam Mattress
Item 650-1501-1633N**

- High density foam base provides a firm, supportive foundation
- Soft foam top layer enhances patient comfort
- Super soft tapered heel section provides added pressure reduction for sensitive heel and ankle
- Soft, breathable top cover moves with the patient to reduce shearing, is fluid resistant and easy to clean
- Designed to meet the requirements of HCPCS code: E1399
- 600 lb maximum weight capacity - EVENLY DISTRIBUTED

Lumex® Deluxe Clinical Care Recliner-Wide (FR587W)

Featuring an increased weight capacity of 450 pounds (Evenly Distributed) and a 24" wide seat. Recliner mechanism features dual gas cylinders that give smooth power assistance for back position adjustment and Trendelenburg activation. Features three positions for patient comfort:

- **First or Upright Position:** Ideal for patient transport and to initiate patient treatment / therapy.
- **Second or First Recline Position:** Ideal for patient treatment, rest, TV and reading.
- **Third or Full Recline Position:** Ideal for patient treatment, rest and reading.

Trendelenburg Position is easily actuated by the caregiver from the back of the recliner via a trigger mechanism. Heat and Massage option available.



FR587W



LiquiCell is a registered trademark of LiquiCell Technologies, Inc.
AEGIS Microbe Shield® is a registered trademark of AEGIS Environments

Lumex® Powered Bariatric Recliner (FR588W)

Featuring a heavy duty welded steel frame and heavy-duty actuators, this recliner can safely position patients of up to 700 pounds EVENLY DISTRIBUTED.

- Features infinite positions with a range for patient comfort via a four-function hand pendant.
- Legrest and backrest can be raised and lowered independently, giving the caregiver complete control over positioning.

**AltaDyne™ XT APM System
Item 774200A**

Superior mattress system that offers both alternating pressure therapy and low air loss. Its patented cell-in-cell construction is designed to prevent bottoming out in the event of a power failure.

Other features include:

- Provides 6 options for mattress pressure, plus an optional static mode
- Dynamic alternating pressure cycles adjustable from 5-20 minutes to provide pressure relief
- 42" wide for larger users
- Includes visible and audible alarms, with Alarm Mute feature
- Highly vapor permeable and fluid resistant cover
- 800 lb maximum weight capacity EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS Code: E0277



- Powered unit is equipped with a rechargeable battery pack for emergency operation during a power failure.
- Quick-release trigger mechanism on push handle allows chair to swiftly move into Trendelenburg Position.
- Drop arm mechanism allows both arms to drop down allowing for easier and safer lateral transfers.



FR588W

Both recliners meet California Technical Bulletin CA TB-133 (Flammability Test)



**Three-in-One Aluminum Commode
Item 2190A**

This Lumex commode has a versatile, economical design that supports up to 400 lb **EVENLY DISTRIBUTED**. The lightweight aluminum frame makes it easy to transport. Other features include:

- Contoured armrests enhance user comfort
- Wide 14" x 16 3/4" deep seat that snaps on easily
- Removable backbar for over the toilet use
- Comes complete with commode pail, cover and splash guard
- Meets the requirements for HCPCS E0168



**Imperial Collection Three-in-One
Drop-Arm Commode
Item 6438A**

A commode with a wide and deep seat surface that is built to last. The commode comes with easy release drop arms for safe lateral transfers. Other features include:

- Extra wide 23" platform seat
- 600 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Cushioned armrests provide added comfort
- Height adjusts in 1" increments from 19 1/2" – 23.5"
- Wide base for added stability
- Complete with commode pail, cover and splash guard
- Meets the requirements for HCPCS E0168

**Bariatric Transfer Bench
Item 7925A**

For the bariatric user in need of an extra wide seating surface. Being adjustable in height, this transfer bench will accommodate a broad range of users. Other features include:

- Drainage holes in seat minimize water build-up
- Height adjustable from 16" to 20" in 1" increments
- 600 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Attractive platinum grey in color



**Bariatric Folding Commode
Item 7109A-2**

Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame.

- Welded steel construction provides strength and stability
- Folds flat for convenient storage and travel (to less than 7")
- Tool-free assembly
- Sturdy plastic seat snaps-on easily
- Comes complete with commode pail, cover and splash shield



**Bariatric Bath Seat
Items 7932A & 7939A**

With a durable blow-molded plastic seat, this bath seat provides added stability in the shower. The anodized aluminum frame is lightweight, durable and rust-resistant. Other features include:

- Height adjustable in 1" increments
- 600 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Available with and without back
- Meets the requirements for HCPCS E0245

- Meets the requirements of HCPCS code: E0168
- Maximum Weight Capacity 650 lb - **EVENLY DISTRIBUTED**

LUMEX®



gf
Proudly Made in the USA



**Lumex® Easy Lift
Item LF1090**

The Easy Lift is designed to meet the needs of the expanding bariatric population up to 600 pounds. Features include:

- Six-point spreader bar with 360° rotation tilts for extra comfort and safety
- Lower base height maximizes stability and fits under most beds (4.7" floor clearance)
- Minimum lift range of 27"; maximum is 71"
- Ergonomically designed hand pendant with pendant hook.
- Warning beep for low battery charge and weight capacity levels

**Imperial Offset Cane
Item 6334A**

With an offset handle, this cane places the user's weight directly over the base for maximum control and balance. Other features include:

- Easy push button height adjustment from 29" to 38"
- Heavy duty steel construction
- 600 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Attractive black finish resists scratching
- Meets requirements for HCPCS E0114



**Mesh Full-Body Sling
Items FM110, FM111, & FM140**

Designed with comfort and safety in mind for the larger user. Quick draining/drying mesh fabric offers full head and neck support. Other features include:

- 450 lb maximum weight capacity
- 600 lb maximum weight capacity **EVENLY DISTRIBUTED** on models FM140 and FMC141
- Available commode sling options (FMC114, FMC115, FMC116 & FMC141)
- Designed to meet the requirements of HCPCS code: E0621

**Bariatric Imperial Steel
Crutches**

Items 3614A & 3615A

A crutch that provides support and comfort for users up to 600 lb maximum weight capacity **EVENLY DISTRIBUTED** capacity. The crutches are made from heavy duty steel yet offer contoured underarm pads for enhanced comfort. Other features include:

- Easy push button height adjustment
- Available in both Adult (5'2" – 5'10") and Tall (5'10" – 6'6") sizes
- Handgrip position and overall height can be easily adjusted



**SURELIFT™ Bariatric Padded Sling
Item UP894**

A sling that provides excellent patient comfort and security. Six point connectivity enhances positioning and safety. Other features include:

- Reinforced back provides additional support
- Height adjustable straps for broad range of users
- Padding on legs and back for comfort
- 600 lb maximum weight capacity **EVENLY DISTRIBUTED**

**Imperial Collection Dual-Release
X-Wide Folding Walker
Item 604070A**

Designed with a wider frame, this walker is ideal for the larger user. The aluminum frame is lightweight and supports up to 600 lb maximum weight capacity **EVENLY DISTRIBUTED**. Other features include:

- Width inside handgrips of 20"
- Dual-release is easy to operate
- Height adjustable from 5'4" – 6'2" for broad range of users
- Also available with wheels (Item 604070W)
- Meets requirements for HCPCS E0148





Walkabout Imperial Four-Wheel Rollator Item RJ4405x

The Walkabout Imperial is ideal for larger users who need a wider ambulatory device. It features a 20" wide padded seat and durable frame that supports up to 450 lb maximum weight capacity **EVENLY DISTRIBUTED**.

Other features include:

- Padded backrest and ergonomic handgrips provide extra comfort
- Easy to operate hand brakes
- Width inside handgrips of 24"
- Robust 8" casters ideal for both indoor and outdoor use
- Comes standard with pouch
- Designed to meet the requirements of HCPCS codes: E0149 & E1056



Bariatric Quad Cane Item 6130A Silver Vein 2/cs

- Offset handle places the user's weight over the base for balance and control
- Low profile, 4-point large base is designed for maximum stability and walking assistance
- Rounded edge base plate helps prevent scraping and bruising of ankles
- Simple height adjustment from 29"- 38" in 1" increments
- 500 lb maximum weight capacity

EVENLY DISTRIBUTED

intended for larger users

- Designed to meet the requirements of HCPCS code: E0105

Everest & Jennings®

Traveler HD Wheelchair Item 3G01xxx series

For users up to 500 lb looking for an affordable solution that is built to last. The Traveler HD comes with double cross-braces, heavy-duty 5/8" diameter axles, and long lasting, double embossed seat and back upholstery.

- Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard)
- Durable carbon steel frame with reinforcement at all critical load bearing points
- Two position axle and caster fork provides 2" seat height adjustment
- 24" x 1" black plastic wheels with black handrims and non-marking polyurethane tires
- Meets requirements for HCPCS K0007



Paramount XD Wheelchair Items 5PX10620, 5PX10630, 5PX10820, & 5PX10830

The Paramount XD is ideal for the bariatric wheelchair user with its offering of 26" and 30" seat widths and a 19 1/2" seat depth. Reinforced welds and gussets on frame make it a durable long lasting solution.

Other features include:

- Black reinforced nylon upholstery
- 8" x 2" solid casters with 24" solid rear wheels
- 650 lb maximum weight capacity **EVENLY DISTRIBUTED**
- Optional swingaway footrest and elevating legrest
- Meets requirements for HCPCS K0007 snap-on tires and integrated handrims
- Chrome steel frame with gusset reinforcements

Bariatric Transport Chair Item EJ777-3

Features:

- Fold down back makes transport and storage easy
 - Black, padded, Nylon upholstery is flame retardant
 - 24" seat width
 - Carry straps on seat allow for ease of lifting and folding
 - Weighs 31 lb without front rigging
 - Weight Capacity 400 lb **EVENLY DISTRIBUTED**
- Meets requirements for HCPCS E1039





**PVC Reclining Shower Chair
Items 89330 & 89440**

Constructed of the highest quality health care grade PVC, to its deluxe elongated extra heavy-duty, contoured open front soft seat provides user with extra comfort. Other features include:

- Ergonomic curved rails provide support without sharp edges with anti-slip handgrips for added safety
- Padded lumbar support and head section – ideal when in use for longer periods
- Dual swing away arms for patient transfer can be secured with red safety straps
- Internal width of 20" (89330) and 26" (89440) with weight capacities of 325 and 600 lb **EVENLY DISTRIBUTED** respectively

PVC Shower Bed Item 8005

A product designed to facilitate patient showering that is both durable and user friendly. This Shower bed comes with a closed cell waterproof foam pad. Other features include:

- Three position elevating headrest
- Side safety railings drop for easy patient transferring
- Full length drain pan catches falling water with hose to drain any residue
- 600 lb maximum weight capacity **EVENLY DISTRIBUTED**

**PVC Shower Chair
Item 89350**

Constructed of the highest quality health care grade PVC to its deluxe elongated extra heavy-duty, contoured open front soft seat provides user with extra comfort. Other features include:

- Quick drying, removable mesh back support is standard in all shower chairs
- Ergonomic curved rails provide support without sharp edges with anti-slip handgrips for added safety
- Optional telescoping footrest feature allows footrest to move in or out with ease. (Item 89351)
- Threaded stem casters provide additional safety
- 700 lb maximum weight capacity **EVENLY DISTRIBUTED**



The Graham-Field® Family of Brands



SIMMONS Clinical Care products manufactured by GF Health Products, Inc. under license from the trademark owner, Dreamwell, Ltd., a subsidiary of Simmons Bedding Company, Atlanta, Georgia. SIMMONS is a trademark of Dreamwell, Ltd. a subsidiary of Simmons Bedding Company. © 2016 Dreamwell, Ltd. All rights reserved.

Information contained herein is subject to change.
The most current and complete product information can be found on our website.



www.grahamfield.com

© 2016, GF Health Products, Inc. All Rights Reserved.
E&J, Everest & Jennings, Graham-Field, and Lumex are trademarks of GF Health Products, Inc.

GF Health Products, Inc. is an ISO 13485 Certified Company.

GF0900037RevD16