

The portable, ultra-compact Nano-Neb Portable Nebulizer is small enough to fit in the palm of the hand, yet powerful enough to deliver fast, effective aerosol treatment. Designed to meet the needs of today's busy lifestyle, this nebulizer is ideal for patients who need mobility and convenience.

- Reliable, Piston-driven Compressor.
- Lightweight – compressor unit weighs only 0.77 lb.
- Mobile - with its small size, quiet operation and optional battery pack, the unit can be used almost anywhere.
- Flexible - can be powered using the AC power cord in a wall outlet, the optional AC car adaptor, or the optional battery pack.
- 3 year limited warranty.
- Designed to meet the requirements of HCPCS Code: E0570 or E0571.

The John Bunn Nano-Neb™ comes complete with:

- 1 reusable nebulizer
- 1 adult mask
- 1 pediatric mask
- 1 PVC tubing
- 1 AC/DC adaptor
- 3 spare filters
- 1 multi-compartment, zippered carry bag

JB0112-065	Nano-Neb Portable Nebulizer	1 ea
JB0112-065B	Nano-Neb with Optional battery + Auto Adaptor	1 ea
JB0112-065BATTERY	Optional battery	1 ea
ADAP3011	Auto Adaptor	1 ea

Optional Accessories (sold separately):

- Portable, rechargeable battery pack
- Auto adaptor



Technical Specifications

Method Of Nebulization	Piston Compressor	Width	3.4"
Power	110-120 VAC, 60Hz	Height	4.9"
Power Consumption	16 watts	Depth	1.9"
Operating Temperature Range	50°F TO 104°F (10°C TO 40°C)	Weight	0.77 lb
Operating Humidity Range	30 to 95% RH	Max Pressure (Compressor)	32.5 psi
Storage Temperature Range	13°F TO 140°F (-25°C TO 60°C)	Particle Size (Nebulizer)	≤ 5 microns
Storage Humidity Range	10 to 95% RH	Capacity (Nebulizer)	7 ml



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com