



The Zenith

The versatile bed that is designed to meet Long Term Care clinical requirements. From a low 7.9 inches to a high position of 30 inches, the Zenith bed facilitates resident safety and ensures caregiver ease of operation. Powered by three ultra-whisper, rapid-moving, DC motors engineered to keep Hi-Lo function in sync. It's the only bed that you'll ever need, and priced so you can afford it!

- Available in 76" and 80" sleep surface lengths.
- 500 lb weight capacity.
- Available with steel-wire grid or pan sleep deck.
- All welds and components are electrostatic powder-coated to provide durable finish.
- Standard with six-function pendant control.
- 115 Volt standard. 230 Volt available.
- UL 60601 and CAN/CSA C22.2 No. 601-1-M90 approved.
- Warranty:
 1. Electronics: 3 Years
 2. All other components: 15 Years


76" Length	80" Length	Sleep Surface	Bed Type	Units
ZG5074	ZG5174	Grid Deck	Full-Electric	1 ea
ZG5014	ZG5114	Pan Deck	Full-Electric	1 ea
ZZ5074	ZZ5174	Grid Deck	Full-Elect. w/ Z-Lok	1 ea
ZZ5014	ZZ5114	Pan Deck	Full-Elect. w/ Z-Lok	1 ea

* Bed, head and foot boards, and siderails sold separately.





Zenith, Full-Electric Bed shown with Z-Lok Pedals, Grid Deck, Custom Windsor Mahogany Center Crown Wood Cap Headboard, Stock Windsor Mahogany Center Crown T-Mold Foot Board with Embedded Staff Control, Head and Foot Half Rails, Assist Bar, and Stabilizer Legs.



Optional Z-Lok Locks at every height

Head and Foot Boards*	
Z46500	Rectangular, Light Oak Finish
Z48500	Rectangular, Medium Oak Finish
Z80500	Rectangular, Dark Walnut Finish
Z46550	Rectangular, Light Oak Finish, Embedded Staff Control**
Z48550	Rectangular, Medium Oak Finish, Embedded Staff Control**
Z80550	Rectangular, Dark Walnut Finish, Embedded Staff Control**
Z29500CCR	Center Crown, Windsor Mahogany Finish
Z29550CCR	Center Crown, Windsor Mahogany Finish, Embedded Staff Control**
Z02500CCR	Center Crown, Figured Mahogany Finish
Z02550CCR	Center Crown, Figured Mahogany Finish, Embedded Staff Control**

* Other finishes available upon request.

** Order in conjunction with ZL74700 Embedded Staff Control.

*** All boards above have T-mold edges.

Side Rails and Assist Bars	
ZA7020	Half Head Rails (Set) or Foot Siderails
AZ7350	Assist Bar (1 Each)

Wall Saver Bumpers	
Z70300	Wall Saver Bumpers, 76" or 80" Bed
554-0766-000	Side Wall Saver Kit

Factory-Installed Options	
ZL71510	Single-Function Lock-Out
ZL71530	Three-Function Lock-Out
ZL74700	Embedded Staff Control
999-0711-005	Grounded Hospital-Grade Plug
554-0748-000	Edema Foot Ratchet

Accessories	
554-0710-000	Stabilizer Legs
554-0719-000	4" Sleep Deck Extension
999-0711-003	Back-Up Battery Pack
999-0711-002	Back-Up Battery Charger
554-0720-000	Trapeze Adapter
554-0143-000	Trapeze Bar (Patient Helper Bar)
554-0190-000	Panel-Mounted I.V. Holder
554-0769-000	Corner Mattress Retainer Kit
554-0743-000	Pendant Holster
554-0721-000	Quick Release Head or Foot Board Mounting Kit
554-0768-000	Traction Adapter Kit

