

WhisperLite[™] II Semi-Electric Homecare Bed



Now available with more durable, high-impact plastic bed ends! Attractive new plastic bed-ends are better able to withstand the storage, handling, cleaning, and use associated with homecare beds.

Introducing the revolutionary WhisperLite $^{\text{\tiny TM}}$ II Semi-Electric Homecare Bed: One of the quietest and lightest beds on the homecare market.

- The revolutionary motor unit and bed design provide the ultimate in quiet, smooth operation.
- Complete welded steel frame construction.
- The head and foot sections can be raised or lowered simultaneously or independently with one hand control.
- During a power failure the motor unit can lower both the head and foot section with just a single standard 9 volt battery.
- Assembles easily with just five components.
- Designed to meet the requirements of applicable HCPCS codes.
- Total bed weight assembled: 154 lb (bed and motor only).
- UL 73 and CAN/CSA 22.2 No. 68-92 approved.
- Available with Fiberboard Decorative Walnut Bed-ends or High-impact Plastic Bed-ends.
- Maximum patient weight 350 lbs.
- 450 lb evenly distributed total weight capacity including patient, mattress, rails, and accessories.
- Bed Height: Min. 15", Max. 24" (at deck level).
- Lift Range: Head 60°, Knee 30°.
- Overall Bed Size: 87" L x 36" W.



Warranty:

- Two year warranty for defects in workmanship and materials of mechanical components, frame, and electronics.
- During the warranty period, defective items will be repaired or replaced at manufacturer's option at no charge.

Bed shipped in three cartons:

- 1. Foot section and motor.
- 2. Head section and Hi/Lo rod.
- 3. Bed end set.

Fiberboard Decorative Walnut Bed-Ends	High Impact Plastic Bed-Ends	Description	HCPCS
WL0208	WL0208PL	Semi-Electric Bed Only	E0295
WL0208-220V	WL0208PL-220V	Semi-Electric Bed Only - 220 V - International Only	
WL0208-PKG	WL0208PL-PKG	Semi-Electric Bed with 1633-Innerspring Mattress and Full Rails	E0260
WL0208-PKGHLF	WL0208PL-PKGHLF	Semi-Electric Bed with 1633-Innerspring Mattress and Half Rails	E0260
WL0208-PKGNGAP	WL0208PL-PKGNGAP	Semi-Electric Bed with 1633-Innerspring Mattress and No-Gap Rail	E0260
WL0208-BFPKG	WL0208PL-BFPKG	Semi-Electric Bed with 1633-Basic Care Foam Mattress and Full Rails	E0260
WL0208-BFPKGHLF	WL0208PL-BFPKGHLF	Semi-Electric Bed with 1633-Basic Care Foam Mattress and Half Rails	E0260
WL0208-BFPKGNGAP	WL0208PL-BFPKGNGAP	Semi-Electric Bed with 1633-Basic Care Foam Mattress and No-Gap Rail	E0260
WL0208-FFPKG	WL0208PL-FFPKG	Semi-Electric Bed with 1633-Extra Firm Foam Mattress and Full Rails	E0260
WL0208-FFPKGHLF	WL0208PL-FFPKGHLF	Semi-Electric Bed with 1633-Extra Firm Foam Mattress and Half Rails	E0260
WL0208-IFPKGFR	WL0208PL-IFPKGFR	Semi-Electric Bed with 1633-Innerspring Mattress with Foam Sleep Surface and Full Rails	E0260
WL0208-IFPKGHLF	WL0208PL-IFPKGHLF	Semi-Electric Bed with 1633-Innerspring Mattress with Foam Sleep Surface and Half Rails	E0260
WL0208-RFPKG	WL0208PL-RFPKG	Semi-Electric Bed with 1633-Reversible Foam Mattress and Full Rails	E0260
WL0208-RFPKGHLF	WL0208PL-RFPKGHLF	Semi-Electric Bed with 1633-Reversible Foam Mattress and Half Rails	E0260
WL0208-XFIPKGFR	WL0208PL-XFIPKGFR	Semi-Electric Bed with 1633-Extra Firm Innerspring Mattress and Full Rails	E0260
WL0208-XFIPKGHLF	WL0208PL-XFIPKGHLF	Semi-Electric Bed with 1633-Extra Firm Innerspring Mattress and Half Rails	E0260
WL0208-RPKG	WL0208PL-RPKG	Semi-Electric Bed with No Mattress and Full Rails	E0261
WL0208-RPKGHLF	WL0208PL-RPKGHLF	Semi-Electric Bed with No Mattress and Half Rails	E0261



For Placing Orders:

1-800-347-5678 Fax: 1-800-726-0601 www.grahamfield.com GF0700058RevB08