



This rigid low bed is designed to be safe and affordable. From a low height of 9.5 inches to a high position of 20.25 inches, the **INTRO** bed utilizes innovative patent-pending lifting technology and quiet, powerful DC motors. With a 450 lb lifting capacity, strong tubular board mounting for maximum push/pull support, and a rigid tubular frame, the Intro offers exceptional durability. Additional standard features include Auto-locking in raised positions for peace of mind, Head & Foot Boards and Patient Assist devices that are removable without tools, and four locking casters. Optional Staff Controls and lockouts also help reduce risk of injury from falls.

**Specifications:**

- Conforms to UL 60601-1 and IEC601-2-38; these are the highest Safety Standards for Medical Beds in Long Term Care, including:
  - No pinch points within 200mm (approx. 8") of the bed perimeter.
  - Rolling clearance of 20mm (approx. 0.8") over thresholds.
  - Spash-proof electrical connections.
  - Load capacity and tip-stability testing.

- 450 lb Dynamic Rated Capacity
- Meets FDA Entrapment Guidelines
- 9.5" Low Height
- 11" of Hi-lo travel
- Rigid Tubular Frame
- Locked when raised and rolls at low height
- Locking Casters all four corners w/360° mobility
- Rigid Tubular Board Mounting
- Head & Foot Boards removable w/o tools
- 85.5" Overall Length with HD/FT Boards
- 39.75" Overall Width w/Rails or Assist devices
- 35" Mattress Deck Width
- Back Mattress Deck articulates to approx. 70 degrees
- Knee Mattress Deck articulates to approx. 30 degrees

76" LENGTH	80" LENGTH	Mattress Deck Surface	Bed Type	Units
NT6074	NT6174	Grid	Full-Electric, 3" Twin Wheel Locking Casters	1 ea
NT6014	NT6114	Pan	Full-Electric 3" Twin Wheel Locking Casters	1 ea
NT6077	NT6177	Grid	Semi-Electric, Fix HT, 3" Twin Wheel Locking Casters	1 ea
NT6017	NT6117	Pan	Semi-Electric, Fix HT, 3" Twin Wheel Locking Casters	1 ea

*Bed, head and foot boards, and rails sold separately.*





*Intro Bed Shown with Stock Windsor Mahogany Center Crown Head and Foot Boards.*

**Options and Accessories:**

Head and Foot Boards*	
<b>Z53500</b>	Rectangular, Solar Oak Finish
<b>Z48500</b>	Rectangular, Medium Oak Finish
<b>Z80500</b>	Rectangular, Montana Walnut Finish
<b>ZL788053</b>	Rectangular, Solar Oak Finish, Embedded Staff Control**
<b>ZL788048</b>	Rectangular, Medium Oak Finish, Embedded Staff Control**
<b>ZL788080</b>	Rectangular, Montana Walnut Finish, Embedded Staff Control**
<b>Z29500CCR</b>	Center Crown, Windsor Mahogany Finish
<b>ZL788029CCR</b>	Center Crown, Windsor Mahogany Finish, Embedded Staff Control**
<b>Z02500CCR</b>	Center Crown, Figured Mahogany Finish
<b>ZL788002CCR</b>	Center Crown, Figured Mahogany Finish, Embedded Staff Control**

Staff Control Electronics Pkg w/o Boards	
<b>ZL788000</b>	For use w/ custom boards

Assist Devices and Assist Bars		
76" LENGTH	80" LENGTH	Descriptions
<b>ZA78400</b>	<b>ZA78400</b>	Half Length Assist Head End (1Pr)
<b>ZA79400</b>	<b>ZA78500</b>	Half Length Assist Foot End (1Pr)
Head End	Foot End	Assist Bars
<b>ZA78600</b>	<b>ZA78650</b>	Assist Bar, 76" & 80" (ea)

Wall Saver Bumpers	
<b>ZA79800</b>	Wall Saver Bumper-Head End

Accessories	
<b>ZA78900</b>	Economy Lockout Box
<b>ZA71900</b>	4" Mattress Deck Extension (supersedes 554-0719-000)
<b>ZA79000</b>	Trapeze Adapter (For use with ZA78100)
<b>ZA78100</b>	Trapeze Bar Offset (1-piece) (must use ZA79000 Adapter)
<b>ZA76800</b>	Traction Adapter Kit (supersedes 554-0768-000)
<b>999-0711-003</b>	Back-Up Battery Pack
<b>999-0711-002</b>	Battery Charger
<b>554-0769-000</b>	Corner Mattress Retainer Kit
<b>554-0190-000</b>	Board Mounted I.V. Sockets

