

## Walkabout Four-Wheel Hemi Rollator

The Lumex Walkabout Four-Wheel Hemi rollator is similar to the Walkabout Lite as it also offers a comfortable, lightweight and affordable solution to mobility. The Walkabout Hemi is ideal for the more petite or shorter individual. Designed with an aluminum frame, this lightweight rollator weighs 14 lbs and supports up to a 300 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.

## Other features include:

- Designed with a lower seat to floor height for shorter users
- Ergonomic hand grips
- Easy to operate locking and quick stop loop brakes (pull up to stop; push down to lock)
- · Padded seat for added comfort
- Removable padded straight backbar
- Adjustable handle height to accommodate broad range of user heights
- 6" wheel for indoor and outdoor use
- Basket included that can be placed underneath seat or in front for carrying of personal items

• Designed to meet the requirements of HCPCS codes: E0143 & E1056

1 ea

| And now wi | th added rear ref | lectors t |
|------------|-------------------|-----------|
| RJ4302B    | Royal Blue        | 1 ea      |
| RJ4302K    | Black             | 1 ea      |

Burgundy

• Lifetime Limited Warranty

## **Specifications**

**RJ4302R** 

| Width Inside Handgrips | 17"    |
|------------------------|--------|
| Overall Width          | 23.5"  |
| Overall Depth          | 21.5"  |
| Seat Width             | 14"    |
| Seat Depth             | 13.5"  |
| Seat Height            |        |
| Caster Diameter        |        |
| Minimum Handle Height  | 32"    |
| Maximum Handle Height  | 36.5"  |
| Weight Capacity        | 300 lb |
| Product Weight         |        |
|                        |        |

Specifications are subject to change without notice







