

Basic Care Foam Mattress**Latex Free**

- Our most economical pressure-reducing mattress offers a high level of therapeutic comfort at an affordable price
- High quality foam base with soft foam comfort layer
- Super soft tapered heel section provides added pressure reduction for sensitive heel and ankle
- Breathable top cover moves with the patient to reduce shearing
- Soft and durable Chemlux II top covers are bacteria resistant, flame retardant, fluid resistant, anti-fungal, anti-static, non-allergenic, have near 0% vapor transmission, and are easy to clean
- Soft and durable Chemsafe bottom covers are bacteria resistant, flame retardant, fluid resistant, anti-static, non-allergenic, stain resistant, and easy to clean
- Meets 16 CFR 1633 requirements
- Designed to meet the requirements of HCPCS code: E0184
- 250 lb maximum weight capacity
- Two year limited warranty
- Latex Free

31580-1633 80" x 36" x 5 ½"**31584-1633** 84" x 36" x 5 ½"