



MATRIX II BED ORDER FORM

999-0807-185
REV E 04/11

Model 709

To order NEW beds, please contact a GF Health Products Sales Representative at 1-800-554-9215

MATRIX II BEDS AVAILABLE			
Order #	Description	Qty Ordered	Order Date
MP Series with Standard Caster Carriage Assembly			
MP7074	76" Grid Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor		
MP7174	80" Grid Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor		
MP7014	76" Pan Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor		
MP7114	80" Pan Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor		
MP Series with Total Lock Caster Carriage Assembly			
MP707404	76" Grid Deck — Full electric with 3 Motors and 4 Total Lock Casters		
MP717404	80" Grid Deck — Full electric with 3 Motors and 4 Total Lock Casters		
MP701404	76" Pan Deck — Full electric with 3 Motors and 4 Total Lock Casters		
MP711404	80" Pan Deck — Full electric with 3 Motors and 4 Total Lock Casters		
MM Series with One Step Locking Caster Carriage Assembly			
MM7074	76" Grid Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor and M-Lok		
MM7174	80" Grid Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor and M-Lok		
MM7014	76" Pan Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor and M-Lok		
MM7114	80" Pan Deck — Full electric with 1 Head, Foot, and Hi/Lo Motor and M-Lok		

NOTE: Footboard with embedded Staff/Nurse Control is optional and is not part of the standard Matrix II Bed.

TABLE OF CONTENTS:

For new or replacement bed service parts please refer to the following pages.

Main Bed Section (Decking, u-levers, radial arm tube, caster carriage, casters, end caps, etc.)	Page 2
Board Components & Boards With Staff Control Assemblies (Boards, Staff Control, Cables, Shroud, etc.)	Page 3
36" Bed End Components (Standard boards, embedded footboards, mounting kits, etc.)	Page 3
Wallsavers and Mattress Retainers	Page 3
Electronic Components (Cables, motors, control boxes, hand controllers/pendants, etc)	Page 4
Economy Lock Out & Battery Backup Options	Page 4
3" Twin Replacement Casters and 4" Total Lock Casters	Page 4
Miscellaneous Hardware and Literature	Page 4
Assist Devices (Half rails, foot mounting bracket, assist bars, trapeze adaptors, trapeze supports, I.V. mount, and traction adaptor)	Page 5
Wide Deck Kits (without boards, with boards and cutouts, and with embedded staff control)	Page 6
Standard 4" Extension Pan Kit and Wide 4" Extension Pan Kit	Page 7
Bed Hardware	Page 8
Control Box and Motor Mounting Hardware	Page 9
Standard Matrix II Bed Diagram	Page 10



MATRIX II BED SERVICE PART

ORDER FORM

999-0807-185
REV E 04/11

Model 709

The following list of service parts for the Matrix II Bed have been provided for your convenience. Please keep this list in a secure location and use as a reference when ordering new or replacement parts. A general bed diagram has been provided on the last page.

To order new or replacement SERVICE PARTS, please contact a Basic American Medical Products Customer Service Representative at 1-800-365-2338 to confirm you are ordering the right part.

MAIN BED SECTIONS		Please see back diagram (page 10) for general part identification.			
Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
DECKING					
Same as Part#	999-0807-900G	Main Frame Weldment (for reference only - not a service part)			
Same as Part#	999-0779-920G	Standard Head/Back Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0779-930G	Standard Foot Grid Deck Weldment (80")	1		
	999-0775-936G	Optional Foot Grid Deck Weldment (76")	(1)		
Same as Part#	999-0779-940G	Standard Knee Grid Deck Weldment (76" & 80")	1		
Same as Part#	999-0796-982G	Optional Head/Back Pan Deck Weldment (76" & 80")	(1)		
Same as Part#	999-0796-983G	Optional Foot Pan Deck Weldment (80")	(1)		
	999-0796-986G	Optional Foot Pan Deck Weldment (76")	(1)		
Same as Part#	999-0796-984G	Optional Knee Pan Deck Weldment (76" & 80")	(1)		
Same as Part#	999-0807-036G	Leveler Strap Bracket	2		
HEAD & FOOT U-LEVERS & DRAG LINK					
Same as Part#	999-0807-913G	Head U-Lever Weldments	1		
Same as Part#	999-0807-903G	Foot U-Lever Weldment	1		
Same as Part#	999-0807-906G	Drag Link Weldment	1		
CASTER CARRIAGE ASSEMBLIES					
Same as Part#	999-0807-955G	Standard Caster Carriage Assembly with two Locking and two Non-Locking Casters	1		
Same as Part#	999-0807-958G	Optional Locking Caster Carriage Assembly with new M-Lok at Foot end and Caster Keepers at Head	1		
Same as Part#	999-0779-954G	Optional Caster Carriage Assembly with four 4" Locking Steel Casters	1		
CASTERS (if ordered separately)					
Same as Part#	100-6330-061	Single 3" Non-Locking Caster-Standard or Lock Carriage	2/4		
Same as Part#	100-6330-062	Single 3" Locking Caster for Standard Carriage	2		
Same as Part#	100-6340-015	Single 4" Total Lock Steel Caster for Total Lock Carriage	4		
CAPPING (if ordered separately)					
Same as Part#	999-0775-001	Half Moon Caps (Head & Foot U-Levers)	12		
Same as Part#	100-4700-015	1.00" Sq. Tube End Cap (Main Frame Crossbraces)	4		
Same as Part#	100-4700-017	1.00" x 2.00" Rect. Tube End Cap (Drag Link & Carriage)	6		

HD/FT BOARDS WITH EMBEDDED STAFF CONTROL ASSEMBLIES

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
ZL788053	-	Solar Oak Hd/Ft. Boards with Staff Control	1/1/1		
ZL788048	-	Med. Oak Hd/Ft. Boards with Staff Control	1/1/1		
ZL788080	-	Montana Walnut Hd/Ft. Boards with Staff Control	1/1/1		
ZL788029CCR	-	Windsor Mahogany Hd/Ft. Boards with Staff Control	1/1/1		
ZL788002CCR	-	Figured Mahogany Hd/Ft. Boards with Staff Control	1/1/1		
ZL788004CCR	-	Wild Cherry Hd/Ft. Boards with Staff Control	1/1/1		

STANDARD BED END COMPONENTS

Please see back diagram (page 10) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
STANDARD HEAD- & FOOTBOARDS WITHOUT CUTOUT:					
Z53500	-	Solar Oak Hd/Ft Boards without staff control cutout	1/1		
Z48500	-	Med Oak Hd/Ft Boards without staff control cutout	1/1		
Z80500	-	Walnut Oak Hd/Ft Boards without staff control cutout	1/1		
Z29500CCR	-	Wind Mah CC Hd/Ft Boards without staff control cutout	1/1		
Z04500CCR	-	Wild Cher CC Hd/Ft Boards without staff control cutout	1/1		
Same as Part#	100-7025-004	1/4-20 UNC x .60 Long Threaded Inserts (pre-inserted into footboard and headboard)	8 (4 per brd)		
STANDARD HEAD- & FOOTBOARDS WITH STAFF CONTROL CUTOUT:					
Z53650	-	Solar Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z48650	-	Med Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z80650	-	Walnut Oak Hd/Ft Boards cutout for emb. staff control	1/1		
Z29650CCR	-	Wind Mah CC Hd/Ft Boards cutout for emb. staff control	1/1		
Z04650CCR	-	Wild Cher CC Hd/Ft Boards cutout for emb. staff control	1/1		
999-0746-888	-	1/4-20 UNC x .60 Long Threaded Inserts (pre-inserted into boards - for replacement order kit of 8)	8 (4 per brd)		
STANDARD MOUNTING KIT 999-0775-000 INCLUDES:					
Same as Part#	999-0746-035G	Angled Mounting Tubes (for mounting boards to frame)	4		
Same as Part#	100-4700-015	1" Square Tube End Caps (pre-inserted in mounting tubes)	4		
Same as Part#	100-5425-094	1/4-20 x 4mm JCB-B Hex Drive bolts	8		
Same as Part#	100-7931-007	5/16 x 1.50" Clevis Pin with hole (mounting tubes to frame)	4		
Same as Part#	100-2001-006	1-9/16 Hair Pin Clip (mounting tubes to frame)	4		
Same as Part#	999-0775-195	Assembly Instruction Sheet	1		

OTHER COMPONENTS

Please see back diagram (page 10) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
-	999-0779-180	Standard Wireform Head-end Wallsaver (Works with 76", 80", 84" & trapeze adaptor bracket positions- also works on Locking Carriage Assembly: 999-0807-958G)	1		
999-0814-000	-	Corner Wire Mattress Retainer Kit (set of two retainers)	1 kit		
Same as Part#	999-0814-001	Single Wire Mattress Retainer	2		
999-0810-000	-	2 Corner Wire Mattress Retainers for Wide Kit Beds	1 kit		

ELECTRICAL COMPONENTS

Please see back diagram (page 10) for general part identification.

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	999-0775-052	Head Actuator/Motor	1		
Same as Part#	999-0778-204	Head Actuator/Motor Cable (8")	1		
Same as Part#	999-0775-053	Foot/Knee Actuator/Motor	1		
Same as Part#	999-0775-201	Foot/Knee Actuator/Motor Cable (60")	1		
Same as Part#	999-0808-950	Hi/Lo Actuator/Motor Assembly	1		
Same as Part#	999-0775-202	Hi/Lo Actuator/Motor Cable (28")	1		
Same as Part#	999-0779-300	J-Box/Control Box	1		
Same as Part#	999-0775-208	J-Box/Control Box 3-Prong Power Cable with Ground	1		
Same as Part#	999-0775-305	Standard Pendant/Hand Controller	1		
Same as Part#	999-0791-000	Pendant Holster	1		
Same as Part#	999-0779-200	Optional Control Box/Staff Control Extension Cable (24")	1		
-	999-0775-901	Optional Staff Control Assembly (with 2 cables, mounting bezel, and control panel)	1 Assm.		
Same as Part#	999-0775-006G	Optional Wide Shroud (covers staff control cables on footboard)	1		
Same as Part#	100-1100-024	4" Cable Tie (staff control cables to footboard mounting tube)	3		
Same as Part#	100-5506-002	#6 x .50" Phillips Truss Head Screws (shroud & control)	6		
ZL788000	-	Electronics Kit for embedded staff control - includes Control Panel, Bezel, Attached Cables, Shroud, Instruction Sheet, attaching screws, and 4" cable tie	1 kit		
Same as Part#	100-1100-018	Cable Ties for all cables to frame	22		

ECONOMY LOCK OUT AND BACKUP BATTERY CHARGER

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
ZA78900	-	Economy Lockout Box (Includes mounting bracket & hardware)	1		
Same as Part#	999-0711-002	Battery Charger/Transformer (order with battery pack)	1		
Same as Part#	999-0711-003	Backup Battery Pack (order with charger/transformer)	1		

REPLACEMENT 3" TWIN CASTERS

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	100-6330-061	Single Non-Locking 3" Twin Casters for Std. Caster Base	2 per base		
Same as Part#	100-6330-062	Single Locking 3" Twin Casters for Std. Caster Base	2 per base		
Same as Part#	100-6330-061	Single Non-Locking 3" Twin Casters for M-Lok	4 per base		

REPLACEMENT HARDWARE FOR TOTAL LOCK 4" CASTER CARRIAGE ASSEMBLY: 999-0779-954G

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	100-6340-015	4" Double Lock Caster with .513" Dia. x 1.625" Stem	4		
Same as Part#	999-0806-150	4" Caster Mounting Bracket	4		
Same as Part#	100-6750-005	1/2-13 High Hex Head Nylon Lock Nut	4		

MISCELLANEOUS HARDWARE & INSTRUCTION LITERATURE

Order #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
Same as Part#	100-8600-015	1.00" Dia. x .25"Thk Duopreen Cushion Disk (head deck)	2		
Same as Part#	100-4550-010	Button Stops (bottom of main frame brackets)	4		
Same as Part#	100-0510-001	Pendant Alligator Clip (for pendant cable)	1		
Same as Part#	999-0807-190	Service/Instruction Manual for Matrix II Bed	1		

ASSIST DEVICES

NOTE: Rail "Kits" include both Quick Release Lanyard and Locking Bolt/Nut mounting hardware

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA78400 (both rails & hardware)	-	Left - Head Half Side Rails (76" & 80") Right - Head Half Side Rails (76" & 80")	1/1/4 for head end		
ZA78500 (includes both rails & two foot bracket kits)	-	Left - Foot Half Side Rails (80") with hardware Right - Foot Half Side Rails (80") with hardware Foot Bracket Kits with hardware	1/1/4 for 80" foot end		
ZA79400 (includes both rails & two foot bracket kits)	-	Left - Foot Half Side Rails (76") with hardware Right - Foot Half Side Rails (76") with hardware Foot Bracket Kits with hardware	1/1/4 for 76" foot end		
	999-0784-907G	Rt. Head (76"/80") or Lft. Foot (80") Half Rail only	1		
	999-0784-908G	Lft. Head (76"/80") or Rt. Foot (80") Half Rail only	1		
	999-0794-907G	76" Bed Left Foot Half Rail Only	1		
	999-0794-908G	76" Bed Right Foot Half Rail Only	1		
-	999-0785-905	Half Rail Foot Bracket Hardware Bag Only (Includes 2 brackets and mounting hardware)	1		
-	999-0775-025G	Foot Bracket only	1		
Assist Bars					
ZA80300 (includes hardware)	-	One Head Assist Bar Kit (Includes Assist Bar and Mounting Hardware)	1/2 for head end		
Same as Part#	999-0803-905	Head Assist Bar Hardware Bag Only (Includes Quick Release Lanyard and Bolt/Nut Hardware)	1 Bag		
ZA80350 (includes 1 foot bracket)	-	One Foot Assist Bar Kit (Includes Assist Bar Kit & Foot Bracket Kit w/mounting hardware)	1/2 for foot end		
Same as Part#	999-0786-995	Assist Bar Foot Bracket Hardware Bag Only (Includes 2 brackets and mounting hardware)	1 Bag		
Same as Part#	999-0803-902G	One Hd. or Ft. Assist Bar Only (without hardware)	1		
Same as Part#	999-0784-100	Quick Release Lanyard Assembly Only	1		
Trapeze Support					
ZA79000	-	Trapeze Support Adaptor Kit (Includes Trapeze Support and attaching hardware)	1 for head end		
Same as Part#	999-0790-901G	Trapeze Support Weldment Only	1		
Same as Part#	999-0790-905	Trapeze Support Hardware Bag Only (Includes mounting hardware and assembly instruction sheet)	1 Bag		
ZA78100	-	Trapeze Bar (Patient Helper) (must use with ZA79000)	1		

ACCESSORIES

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA76800	-	Traction Adaptor Kit	1		
ZA19000	554-0190-000	Brd. Mounted I.V. Sockets (includes bracket, cap & hdwr)	1		

42" WIDE DECK KIT WITHOUT BOARDS

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA811000	-	42" WIDE DECK KIT WITHOUT BOARDS	1		
Kit Includes the Following . . .					
-	999-0811-921G	Left-Hand Head Extension Weldment	1		
-	999-0811-922G	Right-Hand Head Extension Weldment			
-	999-0810-004G	Left-Hand Seat Pan Extension Weldment			
-	999-0810-005G	Right-Hand Seat Pan Extension Weldment			
-	999-0810-002G	Left-Hand Knee Extension Weldment			
-	999-0810-012G	Right-Hand Knee Extension Weldment			
-	999-0811-931G	Left-Hand Foot Extension Weldment			
-	999-0811-932G	Right-Hand Foot Extension Weldment			
-	999-0811-950	Hardware Box	1		
999-0811-950 Hardware Box Includes . . .					
-	999-0032-010	Short Self-Wrenching Nut (Head end pan assemblies to foot end pan assemblies)	2 1 each side		
-	100-5431-002	5/16-18 x .375" Hex Washer Head Machine Screw (Head end pan assemblies to foot end pan assemblies)	2 1 each side		
-	100-7931-008	5/16 x 3.50" Clevis Pin w/Hole (Secures pan sides to deck frame sides (4 RH and 4 LH))	8 4 each side		
-	100-2001-006	1/9/16 Hair Pin Clip (Secures long clevis pins to deck frame)	8 4 each side		
-	999-0811-015	Wide Kit Instruction Sheet	1		

WIDE DECK KITS WITH WIDE BOARDS AND STAFF CONTROL CUTOUTS

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
WIDE KITS FOR MATRIX II BEDS WITH BOARDS & STAFF CONTROL ASSEMBLIES					
ZA811053	-	42" Wide Kit for Rail Mount - Solar Oak Bds w/CUTOUT			
ZA811048	-	42" Wide Kit for Rail Mount - Med Oak Bds w/CUTOUT			
ZA811080	-	42" Wide Kit for Rail Mount - Walnut Oak Bds w/CUTOUT			
ZA811004CCR	-	42" Wide Kit for Rail Mount - CherCC Oak Bds w/CUTOUT			
ZA811029CCR	-	42" Wide Kit for Rail Mount - WinCC Oak Bds w/CUTOUT			

WIDE KITS WITH BOARDS AND EMBEDDED STAFF CONTROL ASSEMBLIES

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
WIDE KITS FOR MATRIX II BEDS WITH BOARDS & STAFF CONTROL ASSEMBLIES					
ZL811053	-	42" Wide Kit for Rail Mount - Solar Oak Bds w/Staff Control			
ZL811048	-	42" Wide Kit for Rail Mount - Med Oak Bds w/Staff Control			
ZL811080	-	42" Wide Kit for Rail Mount - Walnut Oak Bds w/Staff Control			
ZL811004CCR	-	42" Wide Kit for Rail Mount - CherCC Oak Bds w/Staff Control			
ZL811029CCR	-	42" Wide Kit for Rail Mount - WinCC Oak Bds w/Staff Control			

STANDARD 4" EXTENSION PAN KITS

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA71900	-	STANDARD 4" EXTENSION PAN KIT includes:	1		
-	999-0719-001G	Standard 4" Extension Pan - Grey	1		
-	999-0719-995	Hardware Bag	1		
Hardware Bag Includes . . .					
-	100-5231-001	5/16-18 x.75" Carriage Bolt (Attaches pan to foot deck from top)	2 1 each side		
-	100-6931-001	5/16-18 Wing Nut (Secures carriage bolts attaching pan to foot deck)	2 1 each side		
-	100-7931-006	5/16 x 1.75" Clevis Pin (Replaces existing clevis pin for mounting tube positioning)	2		
-	999-0719-115 (p1)	4" Extension Pan Instruction Sheet	1		

WIDE 4" EXTENSION PAN KITS

Purchase #	Part #	Description	Qty Needed On Bed	Qty Ordered	Order Date
ZA821000	-	WIDE DECK 4" EXTENSION PAN KIT	1		
Kit Includes the Following . . .					
-	999-0821-001G	Wide Deck 4" Extension Pan - Grey	1		
-	999-0821-995	Hardware Bag	1		
999-0821-995 HARDWARE BAG INCLUDES . . .					
-	100-5225-009	1/4-18 x.75" Carriage Bolt (Attaches pan to foot deck from top)	2 1 each side		
-	100-8525-001	1/4-18 Spring Lock Washer (Sits on carriage bolt between deck bottom and pan knob)	2 1 each side		
-	999-0810-009	Black Pan Knob (Secures carriage bolts attaching pan to foot deck)	2 1 each side		
-	100-7931-006	5/16 x 1.75" Clevis Pin (Replaces existing clevis pin for mounting tube positioning)	2		
-	999-0719-115 (p2)	Wide Deck 4" Extension Pan Instruction Sheet	1		

BED HARDWARE

Purchase #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
DECK TO FRAME MOUNTING HARDWARE					
Same as Part#	100-5431-002	5/16-18 x .375" Hex Washer Head Machine Screw Secures Head & Knee Decks to Pan Seat, Knee to Foot Deck, and Leveler Straps to Foot Deck Leveler Brackets	8		
Same as Part#	999-0032-010	Short Self-Wrenching Nut Secures Head & Knee Decks to Pan Seat and Knee to Foot Deck	6		
Same as Part#	806-0277-008	Long Self-Wrenching Nut Secures Leveler Straps to Foot Deck Leveler Brackets	2		
Same as Part#	100-7950-004	.50" x 2.00" Clevis Pin with hole Secures Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-7518-001	7/16-.50" Circle Cotter Secures Clevis Pin for Head & Knee Decks to Pan Seat Brackets	2		
Same as Part#	100-6731-006	5/16-18 Thin Nylon Lock Nut Secures Carriage Bolts for Leveler Straps to Frame Leveler Mounting Weldments	2		
Same as Part#	100-5231-005	5/16-18 x 1.75" Carriage Bolt Secures Leveler Straps to Frame Leveler Mounting Weldments	2		
U-LEVERS TO MAIN FRAME AND DRAG LINK MOUNTING HARDWARE					
Same as Part#	999-0779-100	1/2-13 x 2.50" Button Head Shoulder Bolt (T55) Secures Upper U-Levers to Main Frame Brackets	4		
Same as Part#	100-8601-001	.516"ID x 1.25"OD x .062"THK Nylon Washer Between Upper U-Levers to Main Frame Brackets	8		
Same as Part#	100-6750-006	1/2-13 Low Nylon Lock Nut Between Upper U-Levers to Main Frame Brackets	4		
Same as Part#	100-7950-005	.50" x 1.75" Clevis Pin with hole Secures Drag Link to Upper U-Levers	2		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures U-lever and Drag Link Clevis Pins	2		
STANDARD CASTER CARRIAGE MOUNTING HARDWARE					
Same as Part#	100-7950-002	1/2-13NC x 2.50" Button Head Cap Screws (T55) Secures Lower U-Levers to Upright Carriage Mounting Brackets	4		
Same as Part#	100-8601-001	1.25"OD x .516"ID x .062"THK Nylon Washer Between Lower U-Levers and Carriage Mounting Brackets	8		
Same as Part#	100-7500-015	1/2-13 Low Nylon Lock Nut Secures Lower U-Lever Cap Screws to Upright Brackets	4		
Same as Part#	100-6744-001	7/16-14 Full Nylon Lock Nut - Standard Carriage Secures Caster Stems to Caster Carriage Mounting Tubes	4		
OPTIONAL CASTER KEEPER HARDWARE - LOCKING CARRIAGE ONLY					
Same as Part#	999-0809-150	Inner Keeper 2 used on Standard Carriage Caster (head end)	2		
Same as Part#	999-0809-151	Outer Keeper 2 used on Standard Carriage Caster (head end)	2		
Same as Part#	100-5425-059	1/4-20 x 2.50" Hex Head Cap Screw Used to secure Inner Keeper to Outer Keeper	2		
Same as Part#	100-6725-003	1/4-20 Thin Nylon Lock Nut Secures Hex head Cap Screw to Inner Keeper and Outer Keeper	2		
MISCELLANEOUS OTHER MAIN FRAME HARDWARE					
Same as Part#	100-5231-003	5/16-18 X 1.50" Carriage Bolt Used in head end main frame rails to guide mounting tubes	2		
Same as Part#	100-6731-006	5/16-18 Thin Nylon Lock Nut Secures head end carriage bolts to main frame rails	2		

CONTROL BOX AND MOTOR MOUNTING HARDWARE					
Purchase #	Part #	Description	Qty On Bed	Qty Ordered	Order Date
CONTROL BOX MOUNTING HARDWARE					
-	Cannot be ordered separately	Control Box Clip (comes with control box) Secures Control Box to Head Motor Chassis	1		
-	Cannot be ordered separately	Control Box Cord Retainer Clip (comes with control box) Secures Cables to Control Box ports	1		
Same as Part#	100-1100-050	17.5" Nylon Cable Tie Secures Control Box to Head Motor Chassis	1		
Same as Part#	100-5712-002	#12 x .50" Phillips Head Self-Tapping Screw Secures Ground Wire to Main Frame Rail	1		
HEAD MOTOR MOUNTING HARDWARE					
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Head Motor Shaft to Head Deck Bracket	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Head Deck & Motor End Mounting Clevis Pins	2		
FOOT MOTOR MOUNTING HARDWARE					
Same as Part#	999-0775-445	12mm Dia. x 45mm Clevis Pin with hole Secures Foot Motor End to Main Frame Crossbrace Bracket	1		
Same as Part#	999-0808-455	12mm Dia. x 55mm Clevis Pin with hole Secures Foot Motor Shaft to Knee Deck Bracket	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Foot Mounting Clevis Pins	2		
HI/LO MOTOR MOUNTING HARDWARE					
Same as Part#	999-0775-445	12mm Dia. x 45mm Clevis Pin with hole Secures Hi/Lo Motor End to Center Drag Link Bracket	1		
Same as Part#	999-7950-004	.50" x 2.00" Clevis Pin with hole Secures Hi/Lo Motor Assembly Arms to Head-end Drag Link Bracket	1		
Same as Part#	100-7500-015	.50" Reg. .062 Rue Ring Secures Hi/Lo Assembly Control Arms Clevis Pin to Drag Link	2		



STANDARD MATRIX II BED SERVICE PART DIAGRAM

Model 709 Series

