

Aluminum Transport Chair

The Everest & Jennings Transport Chair is designed as an affordable and safe patient transfer device. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

Features:

- Stylish color options of Red, Blue, Silver, Green, Pink and Aqua (Bariatric available in Silver only)
- Fold down back makes transport and storage easy
- Black, padded, Nylon upholstery is flame retardant
- 17" and 19" seat widths (24" for Bariatric model)
- Carry straps on seat allow for ease of lifting and folding
- Footrests are industry compatible and offer easy tool free adjustment (standard models only)
- Full length padded arms provide user comfort (Desk arm for Bariatric model)
- Weighs 20 lb (w/o front rigging) (31 lb for Bariatric model)
- Weight Capacity 250 lb (400 lb for Bariatric model)
- Designed to meet the requirements of Medicare code E1038 (E1039 for Bariatric model)

Warranty:

1 year on frame and crossbraces

3 years on frame and crossbraces for Bariatric model



EJ777-2



6 stylish color options

Item Number	Description	Seat	Seat Depth	Seat Height	Back	Overall Depth	Overall Height	Overall Width	Folded Width	Weight
EJ781-1	Lightweight Aluminum - Red	17"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	22 ¼"	10 ¾"	20 lb
EJ782-1	Lightweight Aluminum - Red	19"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	24 ¼"	10 ¾"	20 lb
EJ783-1	Lightweight Aluminum - Silver	17"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	22 ¼"	10 ¾"	20 lb
EJ784-1	Lightweight Aluminum - Silver	19"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	24 ¼"	10 ¾"	20 lb
EJ787-1	Lightweight Aluminum - Blue	17"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	22 ¼"	10 ¾"	20 lb
EJ788-1	Lightweight Aluminum - Blue	19"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	24 ¼"	10 ¾"	20 lb
EJ789-1	Lightweight Aluminum - Green	17"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	22 ¼"	10 ¾"	20 lb
EJ780-1	Lightweight Aluminum - Green	19"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	24 ¼"	10 ¾"	20 lb
EJ790-1	Lightweight Aluminum - Pink	17"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	22 ¼"	10 ¾"	20 lb
EJ791-1	Lightweight Aluminum - Pink	19"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	24 ¼"	10 ¾"	20 lb
EJ792-1	Lightweight Aluminum - Aqua	17"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	22 ¼"	10 ¾"	20 lb
EJ793-1	Lightweight Aluminum - Aqua	19"	16"	19 ¼"	18 ½"	27 ¼"	39 ¼"	24 ¼"	10 ¾"	20 lb
EJ777-2	Aluminum - Silver (Bariatric)	24"	18"	20"	18"	40"	39 ¼"	30 ¼"	12"	31 lb

NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4". *Without Front Rigging.



GRAHAM-FIELD

www.grahamfield.com

© 2010 GF Health Products, Inc.
Atlanta, Georgia 30360 GF1000033RevB10