

## Exam Lamps Series

- **Maximum of 60 Watt** or 13 Watt CFL equivalent
- 115/120V AC. UL listed 3-wire plug with 8 ft grounded cord
- Chrome-plated pharmacy shade with baked white enamel bowl for maximum illumination
- Solid cast SOK-IT-GUARD™ Safety Lock, slip proof adjustable lock
- Height range: 34 ½" – 53½"
- Overall height: 74" (Base, Rod, Gooseneck, Shade)

*Exam Lamps with 3 wire 8 ft. grounded cord/plug*

<b>1697-1</b>	Chrome-plated Base and Shade	1 ea
<b>1697-1-220</b>	Chrome-plated Base and Shade, 220V	1 ea
<b>1697-1C</b>	With auto clutch-lock adjustment	1 ea
<b>1697-1M</b>	With Mobile Base	1 ea
<b>1697-1B</b>	Black Base and Shade	1 ea

*Exam Lamps with 3 wire 8 ft. grounded cord/hospital grade plug*

<b>1697-2</b>	Chrome-plated Base and Shade	1 ea
<b>1697-2M</b>	With Mobile Base	1 ea

