**Latex Free** 



## **Elite Clinical Care Mattress Series**

- This premium support surface offers outstanding comfort, pressure redistribution, and advanced pressure reduction
- Top layer for head and torso is contour cut softer foam, creating individual cells to enhance comfort, pressure relief, and creating air flow for a cooler sleep surface
- Center layer medium density foam provides additional support for top layer
- Bottom layer high density foam provides support to enhance comfort and extend mattress life
- Sloped heel zone reduces heel pressure
- The heel section uses smaller cells to address the sensitive pressure reduction needs of the heels and ankles
- Optional defined perimeter with side egress cut reduces the need for bed side rails
- Zippered cover
- · Breathable top cover moves with the patient to reduce shearing
- Recovery5-Ultra<sup>™</sup> healthcare fabric with four-way stretch top covers are bacteria resistant, flame retardant, fluid resistant, anti-fungal, anti-static, non-allergenic, have near 0% vapor transmission, and are easy to clean
- Soft and durable Chemsafe<sup>™</sup> bottom covers are bacteria resistant, flame retardant, fluid resistant, anti-static, non-allergenic, stain resistant, and easy to clean
- Meets 16 CFR 1633 requirements
- Designed to meet the requirements of HCPCS code: E0184
- 450 lb maximum weight capacity
- Ten year limited warranty on foam; one year warranty on the cover

 91980H-1633
 80" x 35" x 6"

 91984H-1633
 84" x 35" x 6"

 9198042H-1633
 80" x 42" x 6"

 9198442H-1633
 84" x 42" x 6"



