

Patricia™ Sphygmomanometer

- Black enamel manometer
- Color coordinated black cotton cuff
- Latex two (2) tube bladder
- Standard inflation bulb and end valve
- Pressure measurement from 20 mmHg to 300 mmHg
- Lifetime calibration warranty
- Zippered carrying bag is included

Caution: This product contains natural rubber latex, which may cause allergic reactions.

180	Adult	1 ea
180C	Child	1 ea
180I	Infant	1 ea
180T	Thigh	1 ea
180X	Large Adult	1 ea



<u>Cuff Class</u>	<u>Limb Range (inches)</u>	<u>Bladder (inches)</u>	<u>Cuff Dimensions (inches)*</u>
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

* These are subjective measurements and can vary ± 0.25"

