

Traveler® HD

Everest & Jennings introduces the latest in a long tradition of quality wheelchairs. The Traveler HD provides the ideal solution to the need for a heavy duty frame designed to support up to 500 lb users.

The Traveler HD comes with double crossbraces, heavy-duty 5/8" diameter axles, and a long lasting, embossed back and seat upholstery with a double inner liner for durability.

Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard). Meets Medicare K0007 specifications.

The Traveler HD is manufactured with the same commitment to quality Everest & Jennings is known for.

Standard Features

- Embossed black leatherette upholstery with double inner lining for durability
- Hammer tone paint provides a long lasting and durable finish
- Reinforced frame at all critical load-bearing points
- "Tool Free" adjustable front rigging
- Chart pocket standard on back upholstery
- 24" x 1" black plastic wheels with black handrims and non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swing away footrests
- Durable aluminum black footplates
- Available in Detachable Desk and Full Length Armrests
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum weight capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007



Warranty

5 year warranty on sideframe and crossbraces

Traveler® HD		Dimensions and Weights								
		Seat		Back	Overall				Wt.*	
		Diagram 1	Diagram 2	Diagram 3	Diagram 4	Diagram 5	Diagram 6	Diagram 7	Diagram 8	
3G010320	Detachable Desk Arm, Swingaway Footrest	20"	18"	18" or 20"	16 ½"	34"	36"	28 ½"	12 ½"	57 lbs
3G010330	Detachable Desk Arm, Elevating Legrest	20"	18"	18" or 20"	16 ½"	34"	36"	28 ½"	12 ½"	57 lbs
3G010340	Detachable Full Arm, Swingaway Footrest	20"	18"	18" or 20"	16 ½"	34"	36"	28 ½"	12 ½"	57 lbs
3G010350	Detachable Full Arm, Elevating Legrest	20"	18"	18" or 20"	16 ½"	34"	36"	28 ½"	12 ½"	57 lbs
3G010420	Detachable Desk Arm, Swingaway Footrest	22"	18"	18" or 20"	16 ½"	34"	36"	30 ½"	12 ½"	57 lbs
3G010430	Detachable Desk Arm, Elevating Legrest	22"	18"	18" or 20"	16 ½"	34"	36"	30 ½"	12 ½"	57 lbs
3G010440	Detachable Full Arm, Swingaway Footrest	22"	18"	18" or 20"	16 ½"	34"	36"	30 ½"	12 ½"	57 lbs
3G010450	Detachable Full Arm, Elevating Legrest	22"	18"	18" or 20"	16 ½"	34"	36"	30 ½"	12 ½"	57 lbs
3G010520	Detachable Desk Arm, Swingaway Footrest	24"	18"	18" or 20"	16 ½"	34"	36"	32 ½"	12 ½"	57 lbs
3G010530	Detachable Desk Arm, Elevating Legrest	24"	18"	18" or 20"	16 ½"	34"	36"	32 ½"	12 ½"	57 lbs
3G010540	Detachable Full Arm, Swingaway Footrest	24"	18"	18" or 20"	16 ½"	34"	36"	32 ½"	12 ½"	57 lbs
3G010550	Detachable Full Arm, Elevating Legrest	24"	18"	18" or 20"	16 ½"	34"	36"	32 ½"	12 ½"	57 lbs
3G020500	Fixed Full Arm, Fixed Footrest	24"	18"	18" or 20"	16 ½"	34"	36"	32 ½"	12 ½"	57 lbs

