

Deluxe Toilet Seat Riser With Removable Armrests

- Designed for those who have difficulty sitting down or standing up from the toilet
- Adds 3 ½" to toilet seat height, reducing bending effort
- Securely attaches to toilet bowl using bolts to prevent shifting of the seat
- Designed to work with your existing toilet seat and lid to blend in with existing decor
- Available for round or elongated toilet bowl design
- White, heavy duty plastic construction is strong and durable
- Push buttons allow for easy removal of armrests
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0244
- Lifetime limited warranty

6482R-2 Round toilet bowls design - Retail 2/cs

6482ER-2 Elongated toilet bowl design - Retail 2/cs

<u>Specifications</u>	<u>6482R-2</u>	<u>6482ER-2</u>
Width Between Arms:	19.5"	20.5"
Overall Width:	22.5"	23.5"
Seat Height:	3.5"	3.5"
Seat Depth:	17.25"	19"
Maximum Weight Capacity:	300 lb - EVENLY DISTRIBUTED	

