

## Locking Raised Toilet Seat

- Designed for those who have difficulty sitting down or standing up from the toilet
- Adds 4 ½" to toilet seat height, reducing bending effort
- Front clamping mechanism locks securely to the toilet and prevents shifting of the seat
- Wide, contoured seating surface is comfortable
- White, heavy duty plastic construction is strong and durable
- Available with or without padded armrests
- Available in attractive retail package
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0244
- Lifetime limited warranty



Without Armrests

**6486A** 3/cs  
**6486R** Retail package 3/cs

With Removable Armrests

**6487RA-1** Retail Package 1 ea  
**6487RA** Retail Package 2/cs



**Specifications**

Overall Width ..... 20 ¾"  
 Seat Height ..... 4"  
 Seat Depth ..... 16 ½"  
 Width Between Arms (6487RA, 6487RA-1 only) ..... 18 ½"  
 Maximum Weight Capacity ..... 300 lb EVENLY DISTRIBUTED

