

Padded Transfer Bench

- Provides assistance to those who have difficulty stepping over bathtub walls
- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Side support arm aids in lateral transfer
- Seat height is adjustable in 1" increments
- Reversible to accommodate any bathroom
- Attractive platinum grey color
- Comes assembled with tool-free arm and back attachments
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: EO247
- Lifetime limited warranty

7955A		2/cs
7955A-1		1 ea
7956A	with pail and cover	1 ea
7958A	with Tub Clamp, Swing Arm, Pail and Lid	1 ea



7955A

Specifications

Overall Width	29 ¼"
Overall Depth	16 ¾"
Seat Width	25"
Seat Depth	15 ¾"
Minimum Seat Height	17 ½"
Maximum Seat Height	22"
Clearance Beneath Frame	14" - 18 ½"
Height of Back Support Above Seat	16"
Diameter Commode Opening (7956A & 7958A only)	8"
Maximum Weight Capacity	300 lb EVENLY DISTRIBUTED



7956A



7958A

