

Easy Lift STS

- Heavy-gauge steel construction
- Designed to provide quick, safe assistance to those who have trouble standing
- Unique, ergonomic foot pedal base opening
- Large, easy to grip handles for easy maneuvering
- Lower base height maximizes stability and fits under most beds
- Rubber base coating protects furniture and walls
- Adjustable knee pad unit
- Quick release removable foot platform
- Optional buttock strap for patients requiring additional support
- 24V DC motor allows quiet, smooth operation
- Ergonomically designed hand pendant with pendant hook
- Warning beep for low battery charge and weight capacity levels
- Easy access emergency stop button
- White powder-coated finish
- Accessory storage bag included for slings and charger
- Interchangeable battery works with other Lumex Easy Lifts
- Easy access emergency stop button
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS Code E0635

Warranty:

Frame and spreader bar three-year limited
 Actuator two-year limited
 Batteries one-year limited

LF2020 1 ea

Specifications

| | |
|----------------------------------|-----------------------------|
| Floor clearance..... | 4.5" |
| Spreader bar minimum height..... | 33.0" |
| Spreader bar maximum height..... | 68.0" |
| Overall length..... | 40.0" |
| Minimum base width..... | 25.5" |
| Maximum base width..... | 37.0" |
| Maximum Weight Capacity..... | 400 lb - EVENLY DISTRIBUTED |
| Product weight..... | 117.0 lb |

