

Nursing Home/Home Care Innerspring Mattresses

- Economical, innerspring mattress
- Designed specifically for nursing homes and home care use
- Vented and reversible
- Mattress dimensions: 35" width, 7" height
- Vinyl cover is fluid resistant, anti-microbial, flame resistant, stain/tear resistant, anti-static, and hypo-allergenic
- Easy to maintain, comfortable
- Meets 16 CFR 1633 requirements
- One year limited warranty
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED

| | |
|-------------------------|--------------------------|
| GF1500-175-1633 | 75" Mattress |
| GF1500-180-1633 | 80" Mattress |
| GF1500-184-1633 | 84" Mattress |
| GF41657-180-1633 | 80" Mattress, Extra Firm |

