


### Automatic Overbed Table

- Automatic/infinite stops within height travel range
- Modified “H” base allows use with chair or bed
- Top has plastic-laminated wood finish with vinyl edges and base has grey epoxy finish
- Raised edge style top helps to prevent items from sliding off and catches spills
- Order components separately when requesting the finish to match other furniture in the room
- Supports up to 40 lbs. of evenly distributed weight
- Low profile design allows overbed table to be used with limited clearance beds
- One year warranty



SERIES	Table Height	Table Width	Table Depth	Base Width	Base Depth	Distributed Weight Capacity	Shipping Weight
Low Height	21" – 32"	30"	15"	28.5"	15.25"	40 lbs.	21 lbs.
Standard Height	27.75" - 40"	30"	15"	28.5"	15.25"	40 lbs.	22 lbs.

Tabletop Finish	Low Height Overbed Table		Standard Height Overbed Table	
	Standard Top	Raised Edge	Standard Top	Raised Edge
None, Base only	A793000		A797000	
Solar Oak	A793053	A860053	A797053	A859053
Medium Oak	A793048		A797048	
Montana Walnut	A793080		A797080	
Windsor Mahogany	A793029	A860029	A797029	A859029
Wild Cherry	A793004	A860004	A797004	A859004

\* Other finishes also available.

800.554.9215 Fax 800.524.4113  
View these and more products at [www.grahamfield.com](http://www.grahamfield.com)

Basic American Medical Products is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this literature are based on the latest product information available at the time of printing. The most current product information can be found online at [www.grahamfield.com](http://www.grahamfield.com) © 2013, GF Health Products, Inc.