



# Disposable Blood Lancet

## Disposable Blood Lancet

- Sterilized stainless steel
- Microsharp point provides quick and clean puncture
- Precision engineering assures accuracy to 1/100 mm in thickness
- Superior hardness and sharp point
- Straight side spreads blood easily on glass slide
- Individually wrapped and sealed
- Lancet: 1.5 mm (W) x 3.8 mm (L)
- Made in Japan

3554 200/bx



[www.grahamfield.com](http://www.grahamfield.com)

© 2012 GF Health Products, Inc.  
Atlanta, Georgia 30360 GF1200079RevA12