## **Pediatric Balance Beam Scale**

Mechanical pediatric scales provide precision, accuracy and ensure child safety. Durable steel bases and easy-to-clean trays and seats guarantee long lasting reliability.

- Heavy-duty steel base for improved strength and durability
- Sanitary tray with built-in measuring tape provides safety, reads up to 23" x 1/4". Tray made from easy-to-clean ABS plastic
- · Child seat provides safety and comfort and reduces infant movement
- · Rotating poise bars for dual reading of pounds and kilograms. Poise bar lock-in for error-free operation
- · Factory calibrated to ensure accuracy
- Weight Capacity: 130 LB / 65 Kg
- "LB" measurement marking goes up to 130 LB, and "Kg" measurement marking goes up to 65 Kg

## **2 Year Limited Warranty**

**HOM1524** 1 ea



