

# Slide-W-I-D-E™

## Expandable/Retractable Sleep Deck



**BASIC AMERICAN**  
MEDICAL PRODUCTS



- Expandable/retractable deck quickly and easily adjusts to 35", 39", and 42" with two quick-release ball-detent pins on each side.
- Quick adjustment is ideal for moving through narrow doorways.
- All pins are attached to the bed with lanyards, eliminating lost parts.



- Wide head and footboards available in 39" and 42" widths.
- Rails and assist devices are also quick-release without tools on each side.

**Wider sleep surfaces reduce falls by residents accustomed to sleeping on common consumer products, which range from 39" for a twin to 60" for a queen sized bed.<sup>1</sup>**

<sup>1</sup>Fragala G, Perry B, Fragala M, Examining bed width as a contributor to risk of falls from bed in long term care. *Annals of Long-Term Care: Clinical Care and Aging*, 2012; 20(6):35-38.

# Slide-W-I-D-E™

## Expandable/Retractable Sleep Deck

Proudly Made in America



### Zenith 9000APS

Advanced Positioning System

with 35"- 39"- 42" expandable/retractable deck

APS98674 with pedal lock

APS986741 with pedal lock, backlit handset and underbed light

APS986742 with pedal lock, backlit handset and underbed light, onboard back-up battery

### Zenith 9000

with 35"- 39"- 42" expandable/retractable deck

ZZ98674 with pedal lock

### Zenith 7000APS

Advanced Positioning System

with 35"- 39"- 42" expandable/retractable deck

APS78674 with pedal lock

### Zenith 7000

with 35"- 39"- 42" expandable/retractable deck

ZZ78674 with pedal lock

ZG78674 with locking casters

Patents Pending

The Slide-W-I-D-E and accessories are in full compliance with the FDA guidelines for reducing the risk of bed entrapment, "Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment". Details can be found at: [www.fda.gov](http://www.fda.gov).

#### Highest Safety and Quality Standards Unrestricted UL Listing

Built to conform to UL 60601-1, IEC601-2-38, and CAN/CSA C22.2 No. 601.1 the Zenith Series offers the following, ensuring the utmost resident and staff safety:

- No pinch points within 200mm of bed perimeter
- Bed can be used as a transport device because it passes threshold clearance at ALL heights by 20mm
- Splash-proof electronics
- Load capacity and tip-stability tested and approved



**BASIC AMERICAN**  
MEDICAL PRODUCTS

800.554.9215 Fax 800.524.4113

View these and more products at [www.grahamfield.com](http://www.grahamfield.com)

Basic American Medical Products is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this literature are based on the latest product information available at the time of printing. The most current product information can be found online at [www.grahamfield.com](http://www.grahamfield.com) © 2015 GF Health Products, Inc.

GF1300012RevC15