

FoldAway Onyx Walker

The Lumex FoldAway walker offers convenience and ease of use in this newly featured onyx line walker. Designed with bungee-corded legs, the walker is easy to assemble and allows the walker to fold to half its size for convenient storage and transport.





Bungee corded legs

Walker is half its size when folded!

Features:

- Bungee corded legs legs fold easily and attach to walker frame clips for easy storage
- Dual lever release mechanism requires little dexterity and allows user to fold walker sides independently
- Onyx Finish provides a more upscale look and visual appeal
- Compact Design folds down for easy storage and transport
- **Lightweight** aluminum frame design allows for ease of lifting. Weighs only 6 lbs (2.75KG)
- Height Adjustable simple push button adjusts height in 1" increments
- Maximum Weight Capacity 300 lb EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS: E0135
- Three (3) year limited warranty

Available in Attractive Retail Packaging!



700176C 1 ea

Specifications:

Width Inside Handg	rips16.65"
Width Inside Rear L	egs19.7"
Overall Width	22.8"
Depth Front to Rear	16.5"
Minimum Height	32"
Maximum Height	41"
Approx User Height	5'4" - 6'8"
	300 lb
	EVENLY DISTRIBUTED

