

Features:

- Provides assistance to those who have difficulty stepping over bathtub walls
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Drainage holes in seat minimize water build-up
- Side support arm aids in lateral transfer
- Knock Down design is easy to assemble and saves space for transport and storage
- Seat height is adjustable in 1/2" increments
- Reversible back attachment to accommodate any bathroom
- Attractive platinum grey color
- Maximum Weight Capacity
350 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



Specifications:

Overall Width.....	29 1/2"
Overall Depth.....	17"
Seat Width.....	25 1/2"
Seat Depth.....	16"
Minimum Seat Height.....	19 1/4"
Maximum Seat Height.....	23 1/4"
Clearance Beneath Frame.....	16 1/4" - 20 1/4"
Height of Back Support Above Seat.....	14 1/2"
Maximum Weight Capacity.....	350 lb EVENLY DISTRIBUTED



7927KD-1 1 ea
7927KD-2 2/cs