



Walkabout Four-Wheel Contour Deluxe Rollator

The Lumex Walkabout Four-Wheel Contour Deluxe rollator offers a comfortable, high end and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 16.5 lbs and supports up to a 300 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.



Features:

- Ergonomic hand grips
- Easy to operate locking and “quick stop” loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable padded contoured backbar for extra support
- Adjustable handle height to accommodate broad range of user heights
- Durable 8" wheel for indoor and outdoor use
- Basket that can be placed underneath seat or in front for carrying of personal items
- 300 lb maximum weight capacity - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Includes easy to read, retail-friendly user manual
- Five (5) year limited warranty

And now with added rear reflectors that increase visibility at night!



RJ4805B	Royal Blue	1 ea
RJ4805CH	Champagne	1 ea
RJ4805K	Black	1 ea
RJ4805L	Lavender	1 ea
RJ4805R	Burgundy	1 ea



Color options

**Rollators are not intended to be propelled while seated.
Do NOT use this rollator as a wheelchair or transport device.**

Specifications

Width Inside Handgrips.....	16"
Overall Width.....	25"
Overall Depth.....	27.5"
Seat Width.....	14"
Seat Depth.....	13.5"
Seat Height.....	22.5"
Caster Diameter.....	8"
Minimum Handle Height.....	34"
Maximum Handle Height.....	39"
Weight Capacity.....	300 lb
Product Weight.....	16.5 lbs

