

The Lumex Walkabout Wide Four-Wheel Rollator offers a lightweight and affordable solution with now MORE added comfort. Designed with a wider aluminum frame, this lightweight rollator now offers a width inside the handgrips of **18.5"** providing more room than that of traditional rollators without having to go to a bariatric model.

Weighing only 16 lbs, the rollator also supports up to 350 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport.

And now with added rear reflectors that increase visibility at night!



15.5" WIDE seat for MORE Comfort

Now 18.5" WIDE Between Handgrips



### Features:

- Ergonomic hand grips
- Easy to operate locking and "quick stop" loop brakes (pull up to stop; push down to lock)
- Wide padded seat for added comfort
- Removable padded straight backbar with easy release pull knobs
- Adjustable handle height to accommodate broad range of user heights
- 6" wheel for indoor and outdoor use
- Basket included that can be placed underneath seat or in front for carrying of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Includes easy to read, retail-friendly user manual
- Five (5) year limited warranty

<b>RJ4318AQ</b>	Aqua	1 ea
<b>RJ4318B</b>	Royal Blue	1 ea
<b>RJ4318CH</b>	Champagne	1 ea
<b>RJ4318G</b>	Teal Green	1 ea
<b>RJ4318K</b>	Black	1 ea
<b>RJ4318L</b>	Lavender	1 ea
<b>RJ4318P</b>	Pink	1 ea
<b>RJ4318R</b>	Burgundy	1 ea



Color options

### Specifications

Width Inside Handgrips.....	18.5"
Overall Width.....	27.5"
Overall Depth.....	27.5"
Seat Width.....	15.5"
Seat Depth.....	13"
Seat Height.....	22.5"
Caster Diameter.....	6"
Minimum Handle Height.....	32"
Maximum Handle Height.....	37"
Weight Capacity.....	350 lb
Product Weight.....	16 lbs

**Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.**