

Platinum Collection Bath Seat with Arms - Retail

- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Available with or without backrest
- Available with or without support arms
- Tool Free design
- Designed to meet the requirements of HCPCS code: E0245
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty

7921RA with backrest 3/cs
7931RA without backrest 4/cs

Specifications

Overall Width.....	20"
Overall Depth.....	17"
Seat Width.....	20"
Seat Depth.....	12"
Minimum Seat Height.....	17"
Maximum Seat Height.....	20"
Width Between Arms.....	21"
Maximum Weight Capacity.....	350 lb EVENLY DISTRIBUTED

