

Platinum Collection Bath Seats - Retail

- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Available with or without backrest
- Available with or without support arms
- Tool Free design
- Designed to meet the requirements of HCPCS code: E0245
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty

With backrest

7921R	Standard Grey	3/cs
7921R-1		1 ea
7921RB	Steel Blue	3/cs
7921RB-1		1 ea
7921RC	Chocolate	3/cs
7921RC-1		1 ea
7921RGY	Sea Grey	3/cs
7921RGY-1		1 ea

Without backrest

7931R	Standard Grey	4/cs
7931RB	Steel Blue	4/cs
7931RB-1		1 ea
7931RC	Chocolate	4/cs
7931RC-1		1 ea
7931RGY	Sea Grey	4/cs
7931RGY-1		1 ea

With support arms and backrest

7921RA	Standard Grey	3/cs
---------------	---------------	------

With support arms

7931RA	Standard Grey	4/cs
---------------	---------------	------

Specifications

Overall Width.....	20"
Overall Depth.....	19"
Seat Width.....	20"
Seat Depth.....	12"
Minimum Seat Height.....	14"
Maximum Seat Height.....	22"
Max Weight Capacity.....	350 lb
	EVENLY DISTRIBUTED

** Colored bath seats not available with support arms*



Steel Blue



Sea Grey



Chocolate



Standard Grey