



Platinum Collection Bath Seats

- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- · Available with or without backrest
- Comes assembled with tool-free back attachment
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0245
- Maximum Weight Capacity 350 lb EVENLY DISTRIBUTED
- · Lifetime limited warranty

with backrest

7921A-1 1 ea

7921A-4 4/cs

with backrest and seat cutout

7923A 2/cs

without backrest

7931A-1 1 ea

7931A 4/cs



Specifications

Overall Width	20"
Overall Depth	17"
Seat Width	20"
Seat Depth	12"
Minimum Seat Height	14"
Maximum Seat Height	22"
Maximum Weight Capacity350 lb EVENLY DISTRIBU	TED





