

Maxi-Drain Transfer Bench

- Provides assistance to those who have difficulty stepping over bathtub walls
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Drainage holes in seat minimize water build-up
- · Side support arm aids in lateral transfer
- Seat height is adjustable in ½" increments
- Reversible back attachment to accommodate any bathroom
- Attractive platinum grey color
- Comes assembled with tool-free arm and back attachments
- Maximum Weight Capacity 350 lb EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: EO247
- · Lifetime limited warranty

7927A 2/cs

Specifications





