

Bariatric Easy Lift STS

- Heavy-gauge steel construction
- Designed to provide quick, safe assistance to those who have trouble standing
- Unique, ergonomic foot pedal base opening
- Large, easy to grip handles for easy maneuvering
- Lower base height maximizes stability and fits under most beds
- Rubber base coating protects furniture and walls
- Adjustable knee pad unit
- Quick release removable foot platform
- Optional buttock strap for patients requiring additional support
- 24V DC motor allows quiet, smooth operation
- Ergonomically designed hand pendant with pendant hook
- Warning beep for low battery charge and weight capacity levels
- Easy access emergency stop button
- White powder-coated finish
- Accessory storage bag included for slings and charger
- Interchangeable battery works with other Lumex Easy Lifts
- Easy access emergency stop button
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS Code E0635

Warranty:

- Frame and spreader bar three-year limited
- Actuator two-year limited
- Batteries one-year limited

LF2090 1 ea

Specifications

Floor clearance.....	4.5"
Spreader bar minimum height.....	33.0"
Spreader bar maximum height.....	68.0"
Overall length.....	40.0"
Minimum base width.....	25.5"
Maximum base width.....	37.0"
Maximum weight capacity.....	600 lb
	EVENLY DISTRIBUTED
Product weight.....	122.0 lb

