

Deluxe Sit-to-Stand Padded Slings

- Cushioned support panels under the arm and lumbar areas enhance patient comfort
- The unique cross-over strap and non-slip inner pad keep the sling in place so that it does not ride up
- Optional stand Assist Buttocks Strap is available for initial lift assistance
- 6-month limited warranty
- For use with LF2020 and LF2090 STS lifts

		Torso Girth	Weight Capacity
DSLSA1	Small	26" – 43"	400 lbs
DSLSA2	Medium	29" – 48"	400 lbs
DSLSA3	Large	36" – 59"	400 lbs
DSLSA4	XX-Large-Bariatric	50" – 79"	600 lbs

