

Full-Body Mesh Commode Sling

- Quick draining/drying mesh fabric offers full head and neck support
- For use with 4-point spreader option
- Designed to meet the requirements of HCPCS code: E0621
- 6-month limited warranty



Item	Sling Size	Accent Color	Body Width	Body Length	Commode Opening to Top of Back	Maximum Weight Capacity - Evenly Distributed
FMC114	Medium	Red	42 3/4"	50 3/4"	34 3/4"	450 lb
FMC115	Large	Green	44 3/4"	54 1/4"	38"	450 lb
FMC116	X-Large	Blue	44 3/4"	60 1/4"	43"	450 lb
FMC141	XX-Large	Yellow	42"	51"	41 1/2"	600 lb

