

Silver Collection Low Profile Quad Canes

- Offset handle places the user's weight over the base for maximum balance and control.
- Low profile, 4-point base further enhances stability.
- Rounded edge base plate helps prevent scraping and bruising of ankles.
- Attractive finish resists marring and scratching.
- Easy push-button height adjustment with locking collet for added security.
- Height adjusts from 30" to 39" in 1" increments.
- Tool-free left- to right-hand adjustment is made easily.
- Chrome-plated steel base is strong and durable.
- Choice of standard or Ortho-Ease handgrip and silver or bronze finish.
- Includes laminated color operating instructions.
- 300 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: E0105.
- Limited lifetime warranty.

Standard Grip

Large Base (12" x 8")

6118A	Aluminum	4/cs
6122A	Bronze	4/cs

Small Base (8" x 6")

6119A	Aluminum	4/cs
6123A	Bronze	4/cs

Ortho-Ease® Grip

- Grip angle is designed to adjust to comfortably fit the user's hand and wrist angle and to minimize stress on the joints.
- Ortho-Ease® grip is ideal for arthritics and for those needing a wider grip for greater support.
- Comes complete with convenient wrist strap.

Large Base (12" x 8")

6120A	Aluminum	4/cs
6124A	Bronze	4/cs

Small Base (8" x 6")

6121A	Aluminum	4/cs
6125A	Bronze	4/cs



Ortho-Ease® Grip

Quad Canes: Standard Grip



Quad Canes: Ortho-Ease® Grip

