

Fixed Wheels

- Converts standard walker into wheeled walker
- Wheels allow for walker to roll easily for forward mobility
- The large 5" wheels are ideal for navigating difficult terrain such as pavement, grass or carpet
- Eight adjustment holes to accommodate a broader range of user heights
- Fixed wheels attach to the front of walkers only and are recommended for use with Lumex glide brakes
- Comes complete with glide covers
- Designed to meet the requirements of HCPCS code: E0155
- 300 lb weight capacity
- Designed to fit 1" standard walker tubing
- Limited lifetime warranty

*Wheels are for recommended for use on Lumex Walkers
(We cannot guarantee the performance of these wheels on other industry walkers).*

603630A	3" wheel	1 pr
603650A	5" wheel	1 pr
603750A	5" wheel with leg extension	1 pr
603950A	5" wheel (600 lb weight capacity)	1 pr



603630A



603650A