

Bariatric Trapeze

Helps patient safely change position in bed and transfer from bed to other equipment with minimal assistance.

- Rubber bumper guards protect walls
- Handbar location fully adjustable
- Comes with floor stand
- Meets requirements for HCPCS code E0912
- 14mm wrench required for assembly

2940B 450 lbs assisted weight capacity 1 ea
2960B 600 lbs assisted weight capacity 1 ea

Specifications

	2940B	2960B
Assisted Weight Capacity.....	450 lb	600 lb
Overall Height (including fixed bracket).....	76 7/8"	80 5/8"
Outside Base Width.....	28 7/8"	29 1/4"
Width Between Base Legs	27"	25 1/4"
Base Length	40 1/2"	41"
Weight of Product	72 lbs	90 lbs

