


Neo-Infantometer

- Designed for measurement of premature infants, neonates and children up to approximately 24 months.
- Sliding rule provides easy-to-read, color-coded graduations from 0" – 40" in 1/8" increments and 0 – 100cm in 1cm increments.
- Improved spring hinge automatically secures outer flaps at right angles, assuring accuracy.

1334

1 ea



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com