

Battery-Powered Floor Lifts

Lumex floor lifts are designed to provide safer and more secure transfers of patients by caregivers. Lumex lifts are very sturdy and easy to use.



- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Six point spreader bar with 360° rotation accommodates 2-, 4-, or 6-point slings
- Large, easy to grip handles for easy maneuverability
- Intuitive foot pedal leg opening mechanism frees up hands
- Ergonomic hand pendant with hook
- Emergency stop button and manual lowering mechanism
- Warning beep and visual indicator on control box for low battery charge (24V DC motor)
- Accessory storage bag and how-to-use video included
- Optional digital scale DSC250
- Designed to meet the requirements of HCPCS Code E0635



Frame and spreader bar: Three (3) year limited

Actuator: Two (2) year limited Batteries: One (1) year limited



Model	Floor Clearance	Spreader Bar Minimum Height from Sling Attachment	Spreader Bar Maximum Height from Sling Attachment	Overall Length	Minimum Base Width (outside)	Maximum Base Width (outside)	Safe Working Load (EVENLY DISTRIBUTED)	Product Weight
LF1050	4.7"	27.5"	68.1"	40.2"	26.0"	37.8"	400 lb	114.5 lb
LF1090	4.7"	27.2"	70.9"	50.0"	26.0"	43.0"	600 lb	137.0 lb

