

LUMEX®

A Graham-Field® Family Brand

Adjustable Aluminum Canes



Vinyl Grip

Soft Grip



Features:

- Contoured handle is designed to provide a comfortable, secure hand hold
- Attractive aluminum finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Height adjusts from 30" to 39" in 1" increments
- Available with a choice of grips, either durable vinyl or soft foam, and a choice of colors shown below
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0100
- 250 lb maximum weight capacity - EVENLY DISTRIBUTED
- Three (3) year limited warranty

6220A	Vinyl Grip	Aluminum	6/cs
6221A	Soft Grip	Aluminum	6/cs
6222A	Soft Grip	Bronze	6/cs
6223A	Soft Grip	Black	6/cs
6224A	<u>Soft Grip Combo Pack</u>		6/cs
	2 Aluminum		
	2 Bronze		
	2 Black		

graham
field™

www.grahamfield.com

© 2016 GF Health Products, Inc.
Atlanta, GA 30360 All rights reserved.
GF1500049RevA16