

## Features:

- Offset handle places the user's weight over the base for balance and control
- Low profile, 4-point large base is designed for maximum stability and walking assistance
- Rounded edge base plate helps prevent scraping and bruising of ankles
- Attractive silver vein finish resists scratching
- Simple height adjustment from 29" – 38" in 1" increments with locking mechanism for added security
- Tool-free left to right-hand adjustment is made easily
- Intended for larger users,  
500 lb maximum weight capacity - EVENLY DISTRIBUTED
- Comfortable, contoured PVC handgrip
- Designed to meet the requirements of HCPCS code: E0105
- Three (3) year limited warranty
- Large Base (12" x 8")

6130A Silver Vein 2/cs

