

LUMEX®


A Graham-Field® Family Brand

American Pride Offset Canes

Standard Features:

- Show your American & Military pride spirit with these designs
- Offset handle places the user's weight directly over the base for maximum balance and wrist support
- Soft, foam grip is comfortable for the hands
- Easy push-button height adjustment from 30" – 39" with locking collet for added security
- 250 lb maximum weight capacity - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0100
- Three (3) year limited warranty
- Comes complete with convenient wrist strap

5970CA Camouflage 6/cs
5970FL Flag 6/cs



graham
field™

www.grahamfield.com

© 2017 GF Health Products, Inc.
Atlanta, GA 30360 All rights reserved.
GF1600139RevA17