

Features:

- Arms drop independently below seat level for safe lateral transfers
- Easy-to-release mechanism may be operated by patients lacking manual dexterity
- Arms automatically lock in raised position
- Extra wide platform seat is comfortable for long periods of sitting
- Welded steel construction and extra wide base provide added stability
- Cushioned armrests enhance user comfort
- Height readily adjusts from 19 1/2" - 23", in 1 inch increments
- Removable backbar facilitates use as a bedside commode raised toilet seat or toilet safety frame
- Comes complete with commode pail, cover and splash shield
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0168
- Three (3) year limited warranty



6438A 1 ea

Specifications

Width Between Arms.....	25 1/2"
Overall Width.....	31"
Seat Width.....	23"
Seat Depth.....	18 1/2"
Seat Height Range.....	19 1/2" - 23"
Maximum Weight Capacity.....	600 lb
	EVENLY DISTRIBUTED
Replacement Tip.....	9290