



- Preferred combination of pressure redistribution and support
- Modified convoluted topper enhances immersion characteristics and increases airflow
- Tri-laminate construction with soft heel-zone
- Firm bolster with optional 3" raised bolster
- Fluid-resistant, low friction, anti-microbial top cover with vinyl bottom cover and full zipper
- Designed to meet the requirements of HCPCS Code E0184
- 400 lb (450 lb for 42" wide models) maximum weight capacity, EVENLY DISTRIBUTED
- Meets 16 CFR 1633
- Five (5) year limited warranty on mattress, One (1) year limited warranty on cover

**This product is not made with natural rubber latex*

Dimensions	Flush Bolster	Raised Bolster
76" X 35"	C500B-3576	C500RB-3576
80" X 35"	C500B-3580	C500RB-3580
84" X 35"	C500B-3584	C500RB-3584
76" X 42"	C500B-4276	C500RB-4276
80" X 42"	C500B-4280	C500RB-4280
84" X 42"	C500B-4284	C500RB-4284