

Features:

- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Raised H-frame design facilitates over-toilet positioning
- Contoured, textured handgrips offer a safe, secure hand hold
- Dual release folding mechanism enables user to fold walker sides independently
- Adjustable height fits a broad range of user heights
- Accommodates all leg attachments

- Folded depth 4"
- Designed to meet the requirements of HCPCS code: (without wheels) E0135 (with wheels) E0143
- Includes laminated color operating instructions
- Three (3) year limited warranty
- Models with wheels come with rear glide caps to facilitate use on all surfaces

*Walkers come pre-assembled with wheels in box



Without Wheels

616070A Adult 4/cs

5" Fixed Wheels

616270A-4 Adult 4/cs

616370A-4 Junior 4/cs

3" Fixed Wheels

616470A Adult 4/cs

Specifications

	616070A	616270A-4	616370A-4
Width Inside Handgrips	16 1/2"	16 1/2"	16 1/2"
Width Inside Rear Legs	20 1/8"	20 1/8"	18 5/8"
Overall Width	22 3/4"	26"	25 3/8"
Depth Front to Rear	16 1/2"	16 1/2"	15 1/2"
Minimum Height	32 1/2"	34 1/2"	29 1/2"
Maximum Height	41 1/2"	41 1/2"	34 1/2"
Approx User Height	5'4" - 6'2"	5'4" - 6'2"	4'7" - 5'3"
Weight Capacity - EVENLY DISTRIBUTED	300 lb	300 lb	300 lb