

Everest & Jennings introduces the Traveler® L3 Plus as the lightweight wheelchair addition to the Traveler® line, a brand respected for performance and durability. The Traveler® L3 Plus lightweight frame weighs less than 36 pounds making it the ideal K0003 Medicare coded wheelchair built for short and long term use as well as rental needs.

This lightweight wheelchair comes standard with non-marking polyurethane tires, industry compatible front rigging and adjustable anti-tippers along with a detachable 2" seat extension with hardware. It features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers fixed flip-back desk length arm options for easy lateral movement. The Traveler® L3 Plus is manufactured with the same commitment to quality Everest & Jennings is known for.



2" seat extension with hardware comes standard

### Standard Features

- Weighs less than 36 lb without front rigging
- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Two position axle and caster fork provides 2" seat height adjustment
- Hammertone finish and carbon steel frame provide strength and durability
- Black nylon upholstery with standard chart pocket on back
- Includes fastener tab under seat connecting both seat and back upholstery to eliminate gap
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Fixed flip-back padded desk arms
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003



### Warranty

Five (5) year limited warranty



Chart pocket













Industry compatible front rigging



Anti-tippers come standard

\*This product is not made with natural rubber latex

Traveler® L3 Plus		Dimensions and Weights									
		Seat			Back	Overall				Weight	
		Width	Depth	Height	Height	Depth	Height	Width	Folded Width	Without front rigging	Max Capacity
											
<b>3F012120</b>	Flip-Back Desk Arm, Swingaway detachable footrest	18"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
<b>3F012130</b>	Flip-Back Desk Arm, Elevating legrest	18"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	26"	12 ¼"	34 lb	300 lb
<b>3F012220</b>	Flip-Back Desk Arm, Swingaway detachable footrest	16"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	24"	12 ¼"	34 lb	300 lb
<b>3F012230</b>	Flip-Back Desk Arm, Elevating legrest	16"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	24"	12 ¼"	34 lb	300 lb
<b>3F012320</b>	Flip-Back Desk Arm, Swingaway detachable footrest	20"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb
<b>3F012330</b>	Flip-Back Desk Arm, Elevating legrest	20"	16" + 2" extension	18" or 20"	17 ¼"	29"	35 ½"	28"	12 ¼"	34 lb	300 lb

Traveler® L3 Plus – Common Accessories & Parts		HCPCS Code
<b>90763064A</b>	Heel loop, pair	E0951
<b>90763062A</b>	Wheel lock lever extension, pair	E0961
<b>90763SE16</b>	2" seat extension kit for 16" W	E0981
<b>90763SE18</b>	2" seat extension kit for 18" W	E0981
<b>90763SE20</b>	2" seat extension kit for 20" W	E0981
<b>300-505</b>	Positioning belt 27-47" L (auto-style) push button	E0978
<b>92133040</b>	Positioning belt 36" L (airline style) metal buckle	E0978
<b>90763230A</b>	Anti-tippers adjustable, pair	E0971
<b>907633CSE</b>	Seat rail plug (dust cap) for seat extension, two pairs per bag	
<b>B10320F</b>	O2-IV holder	E2208
<b>B10325</b>	O2 cylinder holder	E2208
<b>EJO2BAG07</b>	O2 holder bag	E2208