


Padded Knock Down Transfer Bench

Designed to provide assistance to those having difficulty stepping over bathtub walls.



Features:

- Comfortable, cushioned seat and backrest are water-tight and easy to clean
- Anodized aluminum frame is lightweight, durable and rust-resistant
- Side support arm aids in lateral transfer
- Seat height is adjustable in 1" increments
- Reversible to accommodate any bathroom
- Knock-down design makes transporting easy
- Quick and easy tool-free assembly
- Maximum Weight Capacity
300 lb - EVENLY DISTRIBUTED



| Item# | Description | UOM |
|----------|---------------------|------|
| 7955KD | Padded Seat | 2/CS |
| 7955KD-1 | Padded Seat | 1 EA |
| 7956KD-1 | Padded Commode Seat | 1 EA |

Specifications:

- Overall Width..... 29"
- Overall Depth..... 23"
- Seat Width..... 24"
- Seat Depth 16"
- Minimum Seat Height..... 18"
- Maximum Seat Height..... 22"
- Diameter Commode Opening (7956KD-1 only)..... 8"
- Maximum Weight Capacity300 lb EVENLY DISTRIBUTED