

EXAM ROOM TABLES



All 4200 Series exam tables have metal pass-through drawers rather than the typical plastic, so items can be accessed on both sides.



The stirrups have a five-position continuous lock, the most in the industry. A procedure tray is independent of the leg rest and may be accessed easily at any time. An all-metal construction for the base table provides longevity and recyclability.

- Three Pass Through Drawers
- Two Storage Drawers
- Pull-Out Leg Rest with Integrated Treatment Tray
- Five Position Continuous Lock Stirrups
- Powder Coat Finish - Foghorn
- Pelvic Tilt

Exam Table with Three Pass-Through Drawers and Two Storage Drawers and Stirrups

4200-XX Series

Height: 32 in / 81 cm

Overall Length: 72 in / 183 cm

Width: 27 in / 69 cm

Foot Step: W 18 in / 46 cm X L 17 in / 43 cm

Paper Roll Holder Widths:
18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Drawer Size: W 18 in / 46 cm X L 16.5 in / 47 cm X D 4 in / 10 cm

Weight Capacity: 500 lbs / 227 kg (EVENLY DISTRIBUTED)

Shipping Weight: 330 lbs / 150 kg

Replacement Top: GF4200RT-XX

Options:

4200-EO-XX	Electrical Outlet
4200-DW-XX	Drawer Warmer
4200-TB133-XX	TB133

- Dual Pedestal Base with DC Electromechanical Drive
- Pull-Out Leg Rest with Integrated Treatment Tray
- Power Assisted Back
- Five Position Continuous Lock Stirrups
- Hand Pendant: Hi-Low, Back, Flat, Return to Chair
- Powder Coat Finish - Foghorn
- Pelvic Tilt
- UL Listed

High/Low Power Table with Power Assisted Back, Hand Control and Stirrups

4600-XX Series

Height: 18.5 in / 81 cm to 35.5 in / 90 cm*

Overall Length: 72 in / 183 cm (Leg Rest Extended)

Width: 27 in / 69 cm

Leg Rest: W 18 in / 46 cm X L 17 in / 43 cm

Paper Roll Holder Widths:
18 in / 46 cm, 19 in / 48 cm, 20 in / 51 cm

Weight Capacity: 700 lbs / 317 kg (EVENLY DISTRIBUTED)

Shipping Weight: 395 lbs / 179 kg

Replacement Top: GF4600RT-XX

Options:

4600-FS-XX	Foot Control
4600-PS-XX	Physician Storage
4600-AA-XX	Assist Arms
4600-TB-XX	TB133
4600-DW-XX	Instrument Warming Kit
4600-FL-XX	All above options except TB133

*+/- .5 in

Standard Color Options

Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Note: Due to the many variations in monitors and browsers, color samples may appear different. Swatch Card is available from your representative.

1.770.368.4700

Information contained herein is subject to change.
The most current and complete product information can be found on our website.
www.hausted.com

Follow us on



© 2018, GF Health Products, Inc. All Rights Reserved.
Hausted is a trademarks of GF Health Products, Inc. GF Health Products, Inc. is an ISO 13485 Certified Company.
GF1800075RevC19