



## Features:

- Offset handle places the user's weight directly over the base for maximum balance and control
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Available in a variety of finishes, grip styles and lengths
- Includes laminated color operating instructions
- Grip angle is designed to adjust to comfortably fit the user's hand and wrist angle and to minimize stress on the joints
- Ortho-Ease® grip is ideal for arthritics and for those needing a wider grip for greater support
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Three (3) year limited warranty
- 250 lb maximum weight capacity - EVENLY DISTRIBUTED

<b>6327</b>	Bronze	Standard Length	( 31"-39")	6/cs
<b>6328</b>	Bronze	Extra Long	( 34"-41")	6/cs

