

**Comfort - Functionality - Durability**



*Multiple Positions for Enhanced Patient Comfort*

**Features:**

- The Lumex FR566DG Deluxe Clinical Care Recliner with Drop Arms features three positions for patient comfort. The fourth position, a caregiver-activated Trendelenburg (or Shock) Position, is easily actuated from either side of the recliner by a conveniently located foot pedal
- Innovative, easy to actuate drop-arm feature. Self-storing drop-arms are simple for the caregiver to activate and enable a safer, easier lateral transfer for the patient and caregiver
- The wide, contoured seat and back, along with an integrated footrest, provide an enhanced level of patient comfort
- Wall saver armrest caps increase durability of the armrests and protect walls from marring
- Wide range of fabrics / colors
- Quick Ship available

- Four Tente<sup>®</sup> swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster
- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Side table surface is smooth and features tapered edges for easy cleaning
- 350 lb maximum weight capacity, evenly distributed

**Options:**

- Heat and Massage
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow



Arm shown in down position



**Dimensions**

Overall Height .....	51.5" / 131 cm
Overall Width .....	33" / 84 cm
Overall Depth.....	38" / 97 cm
Width Between Arms.....	23" / 58 cm
Arm Height From Seat .....	8" / 20 cm
Seat Height .....	22" / 56 cm
Seat Width .....	20.5" / 52 cm
Seat Depth .....	19" / 48 cm
Depth Fully Reclined	
In Trendelenburg.....	76" / 193 cm
Width With Both Side Tables	
in Upright Position .....	53" / 135 cm
Maximum Weight Capacity	
Evenly Distributed .....	350 lb / 158.8 kg